

Double Dropln® Sinks

DOMESTIC U.S. SIZES

Sink No.	Outlet(s)	Wgt (lb)	Dimensions (in)				
			Inside Bowl			Outlet Location	
			Length	Width	Depth	X	Y
D30D	Corner	110	39.0	16.6	10.8	3.5	3.5

All Dimensions are nominal. Cutsheets available upon request.

Round Dropln® Sinks

DOMESTIC U.S. SIZES

Sink No.	Outlet	Wgt (lb)	Dimensions (in)				
			Inside Bowl			Outlet Location	
			Length	Width	Depth	X	Y
DRS12	Center	18	12.0	Round	7.8	Center	

All Dimensions are nominal. Cutsheets available upon request.

Special Order Round, Trifacia & Eliptical Dropln® Sinks

DOMESTIC U.S. SIZES

Sink No.	Outlet	Wgt (lb)	Dimensions (in)				
			Inside Bowl			Outlet Location	
			Length	Width	Depth	X	Y
DRS10	Center	14	10.0	Round	7.8	Center	
DH20C	End	65	29.5	Hexagon	7.0	Center	
D99	End	61	31.5	Eliptical	10.0	11.3	4.5

All Dimensions are nominal. Cutsheets available upon request.

Round Undermount Sinks

DOMESTIC U.S. SIZES

Sink No.	Outlet	Wgt (lb)	Dimensions (in)				
			Inside Bowl			Outlet	
			Length	Width	Depth	X	Y
RS05C	Center	21	11.9	Round	7.5	Center	

All Dimensions are nominal. Cutsheets available upon request.

Special Order Round & Trifacia Undermount Sinks

DOMESTIC U.S. SIZES

Sink No.	Outlet	Wgt (lb)	Dimensions (in)				
			Inside Bowl			Outlet	
			Length	Width	Depth	X	Y
H20C	Center	82	30.0	Hexagon	6.9	Center	

All Dimensions are nominal. Cutsheets available upon request.


End & Wall Sinks With Integral Backsplash

DOMESTIC U.S. SIZES


Sink No.	Outlet	Wgt (lb)	Dimensions (in)			Backsplash Height
			Inside Bowl			
			Length	Width	Depth	
E46	Front	77	22.0	13.6	12.5	3.8
E47	Front	77	22.0	13.6	12.5	5.8
W08C	Center	46	14.6	14.6	11.9	4.3

All Dimensions are nominal. Cutsheets available upon request.


Dimensions Key




Double Dropln® Sink




Dropln® Sink Styles



Dimensions Key



Undermount Sink Styles



End & Wall Sink Style

