




PRODUCT DIMENSIONS



Pegboards

Style No.	Pegs	Dimensions (in)	
		Width	Height
1824	21	18	24
2030	32	20	30
2424	28	24	24
2430	39	24	30
3024	39	30	24
3030	46	30	30
3036	56	30	36
3230	53	32	30
3624	51	36	24
3630	60	36	30
3636	68	36	36
4836	92	48	36

Custom sizes and layouts are available

Pegboard Specs

Pegboard Thickness	1"
Pegboard Colors	All Durcon Colors Available
Peg Lengths	6.5" and 8.0"
Peg Diameter	1/2"
Peg Material	Polypropylene
Peg Colors	Black, White

Please Note: Listed dimensions are typical values only. Variations within normal tolerances are possible. Durcon makes no warranty, express or implied that the results described herein will be obtained under end-use conditions.