



WARRANTY INFORMATION

STAINLESS STEEL DRYING RACKS

Canadian Scientific inc. will replace any drying racks supplied by Canadian Scientific that have defects in workmanship for a period of 2 years from original delivery.

In the event that Canadian Scientific inc. installed the original drying rack, we will remove the defective drying rack and reinstall its replacement.

CLEANING & MAINTENANCE

STAINLESS STEEL DRYING RACKS

Routine cleaning can be accomplished by using warm water and a cloth. Dry with a towel or cloth to prevent water spots. Wipe in the directions of the polish lines. Microfiber cleaning cloths are a good option to use when cleaning stainless steel.

For finger prints, spray with undiluted white vinegar and wipe dry.

For tough staining, use a commercial stainless steel cleaner. Read the directions on the stainless steel cleaner and test in an inconspicuous spot. Be sure to rinse thoroughly and towel dry.

CANADIAN SCIENTIFIC inc.

14 Wilton Road, Guelph, ON N1E 7L5 Canada
Tel. (519) 835-2949

info@canadianscientific.ca info@csilabs.ca
www.canadianscientific.ca www.csilabs.ca