

Laboratory Pegboards


Kemresin Laboratory Pegboards

Kemresin Pegboards are 1" thick and available in the sizes listed below. The number of supplied pegs is listed per size. Kemresin pegboards are supplied exclusively with standard polypropylene pegs measuring ½" dia. x 6⅞" L.

Peg Specification Note: Custom or non-standard peg lengths, styles, or attachment methods are not available and cannot be accommodated through special requests or pricing.


- At the end of the base part number, peg color will be selected by choosing **-PEGWH** = white or **-PEGBK** = black.
- After peg color, Kemresin color will be selected by choosing **BK** = Black, **GR** = Grey, **PT** = Putty, **SL** = Slate.
- Kemresin pegboards can also be ordered with Kemresin Drip Troughs and/or Finished Backs.

Example: PBDR180124-PEGBK BK is a 18" x 24" Black Kemresin Pegboard with Black Pegs



18" High	Length	Qty Pegs
PBDR180124-PEG_ _ _ _	24"L	22
24" High		
PBDR240124-PEG_ _ _ _	24"L	28
PBDR240130-PEG_ _ _ _	30"L	35
PBDR240136-PEG_ _ _ _	36"L	43
30" High		
PBDR300124-PEG_ _ _ _	24"L	39
PBDR300130-PEG_ _ _ _	30"L	49
PBDR300136-PEG_ _ _ _	36"L	60
PBDR300142-PEG_ _ _ _	42"L	70
PBDR300148-PEG_ _ _ _	48"L	81
36" High		
PBDR360124-PEG_ _ _ _	24"L	44
PBDR360130-PEG_ _ _ _	30"L	56
PBDR360136-PEG_ _ _ _	36"L	68
PBDR360142-PEG_ _ _ _	42"L	80
PBDR360148-PEG_ _ _ _	48"L	92

Type 304 Stainless Steel Laboratory Pegboards




Pegboard without Drip Trough
30" High x ¾" Thick x 24" Long

X-020008 with Open Back

X-020010 with Finished Back

Peg Size: ½" diameter x 6" L
Peg Material: Polypropylene
Peg Color: White
Number of Pegs: 32




Pegboard without Drip Trough
30" High x ¾" Thick x 30" Long

X-020009 with Open Back

X-020011 with Finished Back

Peg Size: ½" diameter x 6" L
Peg Material: Polypropylene
Peg Color: White
Number of Pegs: 50




Pegboard with Drip Trough
30" High x ¾" Thick x 24" Long

X-020004 with Open Back

X-020006 with Finished Back

Peg Size: ½" diameter x 6" L
Peg Material: Polypropylene
Peg Color: White
Number of Pegs: 32



Pegboard with Drip Trough
30" High x ¾" Thick x 30" Long

X-020005 with Open Back

X-020007 with Finished Back

Peg Size: ½" diameter x 6" L
Peg Material: Polypropylene
Peg Color: White
Number of Pegs: 50