


# The SOLUTION: Fume Hood (60-100 fpm)

APPLICATIONS: K-12; Higher Education; Industrial; Bio Tech / R&D



- LINER MATERIAL:**
- POLY RESIN (standard)
  - 304 STAINLESS STEEL
  - 316 STAINLESS STEEL
  - PVC
  - SPECIAL

- PRE-PIPING:**
- PRE-PIPING UP
  - PRE-PIPING DOWN

- Optional FH Depth:**
- 39.75" deep w/ 36"d epoxy dished top (add -40"D suffix to part number)

**SPECIAL REQUIREMENTS:** \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

## EH-SOL STANDARD FEATURES

- Interior Baffle Arrangement:
  - Fixed baffles (standard)
- Exhaust Collar:
  - Type 304 s/s - 10" Dia. Round (standard)
  - Other \_\_\_\_\_
- Interior Lighting Fixture:
  - Standard T-8 Fluorescent Light w/ Bulbs (standard)
- Sash Frame & Glass:
  - Frameless w/ Laminated Safety Glass (standard)
- Airfoil:
  - Painted Steel (standard)
- Sash Type:
  - Vertical (standard)
- Pre-Piping / Pre-Wiring (standard)
  - All hoods are pre-wired & pre-piped unless other wise noted
- Sash Drive (standard)
  - Chain & Sprocket Sash Drive
- Sight Line (standard)
  - 36-inch high viewing sightline.

## OPTIONAL ACCESSORIES

- Lattice Assembly:
  - None
  - Aluminum Lattice Assembly
  - Stainless Steel Assembly
- Interior Lighting Fixture:
  - Explosion Proof Light
- Fume Hood Alarm Option:
  - Alarm Tel AFA\_500
  - Alarm Tel AFA\_1000
  - VAV Cutout (Customer provides cutout detail)
- Sash Full-Length Finger Lift:
  - Stainless Steel
- Airfoil:
  - Type 304 Stainless Steel
  - Type 316 Stainless Steel
- Trough:
  - None (standard)
  - Painted Steel
  - Type 304 Stainless Steel
  - Type 316 Stainless Steel
- Ceiling Enclosure:
  - Yes (Height: \_\_\_\_\_)  No
- Finished Back:
  - Yes  No
- Sash Type:
  - Combination
  - Horizontal
  - Interlocking
- Additional Options:
 

|                       |                              |                             |
|-----------------------|------------------------------|-----------------------------|
| Blower                | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Sash Stops            | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Push Button Sash      | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Automatic Sash Return | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Sash Locks            | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Sash Alarm            | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
- Service Fixture Type:
  - Ball valve w/ lever handles
  - Needle valve w/ four-arm handle



**AIR MASTER SYSTEMS CORP**  
 6480 NORTON CENTER DR.  
 MUSKEGON, MI. 49441  
 TEL. (231) 798-1111 FAX (231) 798-4000  
 www.airmastersystems.com

CUSTOMER

Qty: \_\_\_\_\_  
 Color: \_\_\_\_\_