



### **AMS ACID STORAGE CABINET**

The Acid Storage Cabinet is a “cabinet-within-a-cabinet” design with fully welded 18ga. interior and exterior units. The outer cabinet is powdered coated inside and out. The inner unit is coated with a Hi-tech performance polymer alloy coating (Plascoat).

The interior of the cabinet contains one fully adjustable shelf and a removable rear access panel (access panel is furnished on 24” units and larger). These items are also coated with a Hi-tech performance polymer alloy coating (Plascoat).

Exterior depths are as follows: 18”deep cabinets has an interior depth of 15-1/2” and the 22”deep cabinet has an interior depth of 19-1/2”. A 2”x2” 12ga. support angle is shipped with each cabinet for countertop installation. Labels are applied at the factory on the left hand door unless otherwise specified.

Doors are full height with louvers and are available with or without a lock. The back of the cabinet has two venting holes and are plugged with 2” hole plugs. Cabinets are available with a fascia panel and low height doors when specified.

All acid cabinets come with a removable moulded polypropylene spill tray for easy/safe cleaning and disposal of spills.

Cabinets are guaranteed against manufacturing defects and corrosion (under normal operating procedures) for 5 years.

