


## Pegboards

DOMESTIC U.S. SIZES

Style No.	Pegs	Dimensions (in)	
		Width	Height
1824	21	18	24
1830	25	18	30
2030	32	20	30
2424	28	24	24
2430	39	24	30
2436	36	24	36
3024	39	30	24
3030	46	30	30
3036	56	30	36
3230	53	32	30
3630	60	36	30
3636	68	36	36
4230	67	42	30
4830	81	48	30
4836	104	48	36

Custom sizes and layouts are available



## Pegboard Specs

DOMESTIC U.S. SIZES

Pegboard Thickness	1"
Pegboard Colors	All Durcon Colors Available
Peg Lengths	6.5" and 8.0"
Peg Diameter	1/2"
Peg Material	Polypropylene
Peg Colors	Black, White

Epoxy resin Pegboards are supplied with oval polypropylene pegs angled for effective support and drainage of lab glassware. As an option, 2" [51mm] wide stainless steel drip trays are available for most pegboard sizes. PVC tubing is included.

