

CSi Canadian Scientific

The Future of Laboratory Design

Stainless Steel Laboratory Casework

PRODUCT CATALOGUE



INTER DYNE SYSTEMS






Since 1985, Canadian Scientific has been involved in the design, supply, and installation of laboratories across Canada. We have supplied countless laboratories for agencies of Canada's Federal and Provincial Governments, Canadian universities, school boards, healthcare facilities as well as laboratories for some of North America's most recognized commercial and industrial companies.



Contents

- Cabinet Specifications.....4
- Ordering Guide.....5
- Options Menu.....6
- Base Cabinets — Standing Height 7
 - Open Cabinets.....7
 - Door Cabinets 7
 - Drawer Cabinets.....8
 - Door & Drawer Cabinets.....8
 - Corner Cabinets.....9
 - Sink & Fume Hood Cabinets.....10
- Mobile Cabinets.....12
- Hanging Cabinets.....12
- Tall Cabinets 13
- Wall Cabinets - 30” Height.....14
- Filler Panels.....15
- Knee Space Rails and Panels.....17
- Accessories 18
- Countertop Wall Bracket.....19
- Wall Mounted Sink Enclosure.....19
- Utility Tables.....20
- Table Tops.....21
- Scullery Sinks.....22
- Canopy Hoods.....23
- Umbilicals 23
- Counter Tops 24
- Free Standing Shelves.....25
- Wall Shelving 26
- Seismic Lips.....28
- Acrylic.....29
- Drip Trough and Accessories.....30
- Stainless Steel Pegboards.....31

Cabinet Specifications



Removable backs

Full extension, ball bearing,
plated drawer slides

Doors and drawer heads are "double
pan" construction and sound deadened

Adjustable doors with full
access opening

Full depth adjustable shelf
Adjustable on 2" centers

Stainless steel, 5 knuckle
institutional type hinge — 13 ga.

All surfaces shall be Stainless Steel. Specific test data is available upon request. Type 304.

Gauges of base cabinet components

Aprons	18 ga.
Back Panels	18 ga.
Bottom Panels	18 ga.
Drawer/Door Outer Panel	20 ga.
Door Inner Panel	20 ga.
Drawer Bodies	20 ga.
Legs - 2"sq. tube	16 ga.
Shelves	18 ga.
Side Panels	18 ga.
Table Frames	18 ga.
Security Panels	20 ga.

For complete specifications, request our full line of product specifications. Subject to change

Ordering Guide

S	TYPE	HEIGHT	WIDTH	STYLE	OPTIONS
Cabinet Type					
B	Base (Inset)				
F	Floor (Tall)				
HB	Hanging Base				
MB	Mobile Base				
W	Wall				
Base Cabinet Height					
29"	Sitting Height				
32"	ADA Height				
35"	Standing Height				
Wall Cabinet Height					
24" 30" 36" 48"					
Cabinet Width					
12" 15" 18" 21" 24" 30" 36" 42" 48"					
Cabinet Style					
Reference CSi catalog for door and drawer configurations					
Option Suffix (for base cabinet)					
18	Depth (contact CSi for approval)		PS	1 Pull Out Shelf	
LK	Locks		2PS	2 Pull Out Shelves	
LTK	Toe Kick On Left Side Of Cabinet		3PS	3 Pull Out Shelves	
RTK	Toe Kick On Right Side Of Cabinet		VD	Vertical Divider	
SS	Stainless Steel, Type 304		FB	Finished Backs	
SS316	Stainless Steel, Type 316		EC	End Covers	
Option Suffix (for wall cabinet)					
16	16" Depth		TG	Tempered Glass	
LK	Locks		LG	Laminate (Safety) Glass	
SS	Stainless Steel, Type 304				
SS316	Stainless Steel, Type 316				
Option Suffix (for mobile cabinet)					
18	Depth (contact CSi for approval)		2"	Casters (in lieu of standard 3")	
LK	Locks				
SS	Stainless Steel, Type 304				
SS316	Stainless Steel, Type 316				
Option Suffix (for tall cabinet)					
LK	Locks		LG	Laminate (Safety) Glass	
SS	Stainless Steel, Type 304		PS	1 Pull Out Shelves	
SS316	Stainless Steel, Type 316				
TG	Tempered Glass				

Options Menu

DOOR & DRAWER FRONTS



Inset Steel
(Standard Default)

MAGNETIC CATCH



Magnetic Catch
(Standard Default)

DRAWER SLIDES



100 lb. Full Extension (Standard Default)

PULLS



Stainless Steel
Wire Pulls

HINGES



Inset Steel
Stainless Steel
5-Knuckle

GLASS FRAMED DOORS



- 3/16" Clear Float Glass
(Standard Default)
- 1.4" Laminated Safety Glass
- 1.4" Tempered Safety Glass

Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 22.375" deep. Units are available in various widths as shown. Door units feature a single shelf as standard. All open units have 2 shelves. Additional shelves are available. Specify left or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".



OPEN CABINETS



Width	Item Number
12"	SB3512-000-SS
15"	SB3515-000-SS
18"	SB3518-000-SS
21"	SB3521-000-SS
24"	SB3524-000-SS

Up to 24" wide



30" and wider

Width	Item Number
30"	SB3530-000-SS
36"	SB3536-000-SS
42"	SB3542-000-SS
48"	SB3548-000-SS

DOOR CABINETS



Example:
SB3515-100R

Width	Left & Right Hinged
12"	SB3512-100L-SS
15"	SB3515-100L-SS
18"	SB3518-100L-SS
21"	SB3521-100L-SS
24"	SB3524-100L-SS
12"	SB3512-100R-SS
15"	SB3515-100R-SS
18"	SB3518-100R-SS
21"	SB3521-100R-SS
24"	SB3524-100R-SS

Note: L(left) or R(right) refers to door hinge side.



30" and wider

Width	Item Number
30"	SB3530-200-SS
36"	SB3536-200-SS
42"	SB3542-200-SS
48"	SB3548-200-SS

Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Door units feature a single shelf as standard. Additional shelves are available. Specify left or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".



DRAWER CABINETS



Width	Item Number
12"	SB3512-040-SS
15"	SB3515-040-SS
18"	SB3518-040-SS
21"	SB3521-040-SS
24"	SB3524-040-SS

Up to 24" wide



Width	Item Number
30"	SB3530-040-SS
36"	SB3536-040-SS
42"	SB3542-040-SS
48"	SB3548-040-SS

30" and wider

DOOR & DRAWER CABINETS



Right-hinged door shown

Width	Left & Right Hinged
12"	SB3512-110L-SS
15"	SB3515-110L-SS
18"	SB3518-110L-SS
21"	SB3521-110L-SS
24"	SB3524-110L-SS
12"	SB3512-110R-SS
15"	SB3515-110R-SS
18"	SB3518-110R-SS
21"	SB3521-110R-SS
24"	SB3524-110R-SS



30" and wider

Width	Item Number
30"	SB3530-220-SS
36"	SB3536-220-SS
42"	SB3542-220-SS
48"	SB3548-220-SS

Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Door units feature a single shelf as standard. Additional shelves are available. Specify left or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".



height

drawer

Sink

Right

Example: SC3514-110-SVR-SS

width

door

half drawer

Louver

CORNER CABINETS



Right-hinged door shown

Left Hinged	Right Hinged
SC3514-100L-SS	SC3514-100R-SS


Lazy Susan Add-On Sufflix "LS"



Left Hinged	Right Hinged
SC3514-100-SL-SS	SC3514-100-SR-SS

Left Hinged with Louvers	Right Hinged with Louvers
SC3514-100-SVL-SS	SC3514-100-SVR-SS

Lazy Susan Add-On Sufflix "LS"



Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Door units feature a single shelf as standard. Additional shelves are available. Specify left or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".



Example: **SSB3552-200-SI-SS**

height drawer Sink

width door half Island

 drawer

SINK & FUME HOOD CABINETS



Right-hinged door shown

Width	Left & Right Hinged
12"	SB3512-100-SL-SS
15"	SB3515-100-SL-SS
18"	SB3518-100-SL-SS
21"	SB3521-100-SL-SS
24"	SB3524-100-SL-SS
12"	SB3512-100-SR-SS
15"	SB3515-100-SR-SS
18"	SB3518-100-SR-SS
21"	SB3521-100-SR-SS
24"	SB3524-100-SR-SS



30" and wider

Width	Item Number
30"	SB3530-200-S-SS
36"	SB3536-200-S-SS
42"	SB3542-200-S-SS
48"	SB3548-200-S-SS



Width	Item Number
52"	SB3552-200-S-SS
58"	SB3558-200-S-SS



Width	Item Number
52"	SB3552-200-SI-SS
58"	SB3558-200-SI-SS

Base Cabinets — ADA Height

Standing height base cabinets are 31.75" high and 22.375" deep. Units are available in various widths as shown.

SINK AND FUME HOODS



Width	Item Number
36"	SADAKSC36-SS
42"	SADAKSC42-SS
48"	SADAKSC48-SS



Width	Item Number
30"	SADAKSP30-SS
36"	SADAKSP36-SS
42"	SADAKSP42-SS
48"	SADAKSP48-SS

Mobile Cabinets

All cabinets are available in a mobile version including low profile (27.18" high), sitting height (28.75" high), ADA height (31.75" high), and standing height (34.75" high). They come standard with a finished back and optional finished top. To specify a mobile cabinet use the prefix "M" to the base cabinet model number (example: SMB2724-110). They come standard with a 3" caster. 2" casters are available by adding the suffix "-2" to the part number. Cabinet depth with finished back is 22.812". All mobile cabinets with drawers come standard with a tip resistant counter weight.

Tops: Please add the suffix "T" for steel, "E" for epoxy, "P" for phenolic. All tops are 3/4" thick and with a 7/8" overhang on front only. (Example: SMB2724-110L-T-)




Hanging Cabinets

All cabinets are available in a hanging version including low profile (23.18" high), sitting height (24.875" high), ADA height (27.875" high), and standing height (30.875" high). To specify a hanging cabinet use the prefix "H" to the base cabinet model number (example: SHB2324-110). Finished backs are optional and can be specified with the suffix "-FB". Low profile cabinets should be used for adjustable height work benches or table frames on casters to allow the caster to swivel if the adjustable leg is positioned in the lowest setting.



Tall Cabinets

Tall cabinets are 84" high and available in 16", 18", 21" and 24" deep by using the two digit suffix. Suffix is not needed when selecting 16" deep. Units are available in various widths as shown. Door units feature five shelves as standard. Additional shelves are available. Specify left or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required use the suffix "LK". Wardrobe cabinets are available by adding the suffixes.



OPEN CABINETS



Width	Item Number
18"	SF8418-000-SS
24"	SF8424-000-SS

Up to 24" wide



Width	Item Number
30"	SF8430-000-SS
36"	SF8436-000-SS
42"	SF8442-000-SS
48"	SF8448-000-SS

30" and wider

DOOR CABINETS



Width	Left & Right Hinged
18"	SF8418-1SHL-SS
24"	SF8424-1SHL-SS
18"	SF8418-1SHR-SS
24"	SF8424-1SHR-SS

Right-hinged door shown



Width	Item Number
30"	SF8430-2SH-SS
36"	SF8436-2SH-SS
42"	SF8442-2SH-SS
48"	SF8448-2SH-SS



Width	Left & Right Hinged
18"	SF8418-1GHL-SS
24"	SF8424-1GHL-SS
18"	SF8418-1GHR-SS
24"	SF8424-1GHR-SS



Width	Item Number
30"	SF8430-2GH-SS
36"	SF8436-2GH-SS
42"	SF8442-2GH-SS
48"	SF8448-2GH-SS

Wall Cabinets — 30” Height

Wall cabinets are available in 24”, 30”, 36” and 48” heights and come standard with a 13” depth. For 16” deep cabinets, use the suffix “16”. Units are available in various widths as shown. They feature two shelves as standard for 30” units. Specify left- or right-hinged doors using the suffix “L” or “R” after the two digit depth number. If locks are required, use the suffix “LK”.



OPEN CABINETS



Up to 24” wide

Width	Item Number
18”	SW3018-000-SS
24”	SW3024-000-SS



30” and wider

Width	Item Number
30”	SW3030-000-SS
36”	SW3036-000-SS
42”	SW3042-000-SS
48”	SW3048-000-SS

STEEL-HINGED DOOR CABINETS



Right-hinged door shown

Width	Left & Right Hinged
18”	SW3018-1SHL-SS
24”	SW3024-1SHL-SS
18”	SW3018-1SHR-SS
24”	SW3024-1SHR-SS



30” and wider

Width	Item Number
30”	SW3030-2SH-SS
36”	SW3036-2SH-SS
42”	SW3042-2SH-SS
48”	SW3048-2SH-SS

FRAMED GLASS HINGED DOOR CABINETS



Right-hinged door shown

Width	Left & Right Hinged
18”	SW3018-1GHL-SS
24”	SW3024-1GHL-SS
18”	SW3018-1GHR-SS
24”	SW3024-1GHR-SS




Width	Item Number
30”	SW3030-2GH-SS
36”	SW3036-2GH-SS
42”	SW3042-2GH-SS
48”	SW3048-2GH-SS

OPEN CABINETS



Right-hinged door shown

Lazy Susan Add-On Suffix “LS”



Item Number	Left Hinged	Right Hinged	Left Hinged	Right Hinged
SWC3018-000-SS	SWC3018-1SHL-SS	SWC3018-1SHR-SS	SWC3018-1GHL-SS	SWC3018-1SHR-SS

Filler Panels

Base cabinet filler panels are available in 29", 32" and 35" high. Floor case filler panels are available in 84" high. They are all available in any width required by entering the dimension in lieu of the "XX" in the part number. **Note:** AMS base cabinets are 22.375" deep. (Example: SLF35-6.625)

C-FILLER



Item Number	Item Number	Item Number
SCF29-XX-SS	SCF32-XX-SS	SCF35-XX-SS

C-FRONT FILLER



Item Number	Item Number	Item Number
SFFC29-XX-SS	SFFC32-XX-SS	SFFC35-XX-SS

L-FILLER



Item Number	Item Number	Item Number
SLF29-XX-SS	SLF32-XX-SS	SLF35-XX-SS

L-FRONT FILLER



Item Number	Item Number	Item Number
SFFL29-XX-SS	SFFL32-XX-SS	SFFL35-XX-SS

Filler Panels

Base cabinet filler panels are available in 29", 32" and 35" high. Floor case filler panels are available in 84" high. They are all available in any width required by entering the dimension in lieu of the "XX" in the part number. **Note:** AMS base cabinets are 22.375" deep. (Example: SLF35-6.625)

CORNER FILLER



Item Number	Item Number	Item Number
SICF29-XX-SS	SICF32-XX-SS	SICF35-XX-SS

FLOOR CASE FILLER



Item Number
SFFL84-XX-SS

TOE SPACE FILLER



Item Number
STSF-XX-SS

WALL CASE FRONT FILLER



Item Number	Item Number	Item Number	Item Number
SWCFF24-XX-SS	SWCFF30-XX-SS	SWCFF36-XX-SS	SWCFF48-XX-SS

Wall cabinet front filler panels are available in 24", 30", 36" and 48" high. They are all available in any width required by entering the dimension in lieu of the "XX" in the part number. To specify the depth use the suffix "-13" or "-16".

WALL CASE CORNER FILLER



Item Number	Item Number	Item Number	Item Number
SWCCF24-SS	SWCCF30-SS	SWCCF36-SS	SWCCF48-SS

Knee Space Rails And Panels

KNEE SPACE RAILS



Item Number
SKSR0418-SS
SKSR0424-SS
SKSR0430-SS
SKSR0436-SS
SKSR0442-SS
SKSR0448-SS
SKSR0452-SS
SKSR0454-SS
SKSR0460-SS
SKSR0466-SS
SKSR0472-SS
SKSR0484-SS
SKSR0496-SS

KNEE SPACE PANELS



The suffix "-2P" represents two piece panels.



The suffix "-3P" represents three piece panels.

Item Number	Item Number	Item Number
SKSP2424-SS	SKSP2924-SS	SKSP3524-SS
SKSP2430-SS	SKSP2930-SS	SKSP3530-SS
SKSP2436-SS	SKSP2936-SS	SKSP3536-SS
SKSP2442-SS	SKSP2942-SS	SKSP3542-SS
SKSP2448-SS	SKSP2948-SS	SKSP3548-SS
SKSP2454-SS	SKSP2954-SS	SKSP3554-SS
SKSP2460-SS	SKSP2960-SS	SKSP3560-SS
SKSP2466-SS	SKSP2966-SS	SKSP3566-SS
SKSP2472-SS	SKSP2972-SS	SKSP3572-SS

Item Number	Item Number	Item Number
SKSP2424-2P-SS	SKSP2924-2P-SS	SKSP3524-2P-SS
SKSP2430-2P-SS	SKSP2930-2P-SS	SKSP3530-2P-SS
SKSP2436-2P-SS	SKSP2936-2P-SS	SKSP35236-2P-SS
SKSP2442-2P-SS	SKSP2942-2P-SS	SKSP3542-2P-SS
SKSP2448-2P-SS	SKSP2948-2P-SS	SKSP3548-2P-SS
SKSP2454-2P-SS	SKSP2954-2P-SS	SKSP3554-2P-SS
SKSP2460-2P-SS	SKSP2960-2P-SS	SKSP3560-2P-SS
SKSP2466-2P-SS	SKSP2966-2P-SS	SKSP3566-2P-SS
SKSP2472-2P-SS	SKSP2972-2P-SS	SKSP3572-2P-SS

Item Number	Item Number	Item Number
SKSP2424-3P-SS	SKSP2924-3P-SS	SKSP3524-3P-SS
SKSP2430-3P-SS	SKSP2930-3P-SS	SKSP3530-3P-SS
SKSP2436-3P-SS	SKSP2936-3P-SS	SKSP3536-3P-SS
SKSP2442-3P-SS	SKSP2942-3P-SS	SKSP3542-3P-SS
SKSP2448-3P-SS	SKSP2948-3P-SS	SKSP3548-3P-SS
SKSP2454-3P-SS	SKSP2954-3P-SS	SKSP3554-3P-SS
SKSP2460-3P-SS	SKSP2960-3P-SS	SKSP3560-3P-SS
SKSP2466-3P-SS	SKSP2966-3P-SS	SKSP3566-3P-SS
SKSP2472-3P-SS	SKSP2972-3P-SS	SKSP3572-3P-SS

Cabinet Finished Backs

Cabinet finished backs come standard with a toe kick. Finished backs can be applied to a multitude of cabinets. The example above represents the layout for the part number. The height, width, and depth must match the cabinet that the finished back is to be applied too. If you do not want a toe kick on your unit be sure to put the suffix “NTK” at the end of the part number.

Cabinet finished back - toe kick



Cabinet finished back - no toe kick



Example: CFB2930-18-NTK

height

depth

width

no toe kick
(toe kicks come standard on units)

Accessories

CORNER SLOPED TOP



Item Number
SST-Corner13-SS
SST-Corner16-SS
SST-Corner21 -SS

SLOPED TOP WITH END CAPS




Item Number	Item Number	Item Number
SST13-18-SS	SST16-18-SS	SST21-18-SS
SST13-24-SS	SST16-24-SS	SST21-24-SS
SST13-30-SS	SST16-30-SS	SST21-30-SS
SST13-36-SS	SST16-36-SS	SST21-36-SS
SST13-42-SS	SST16-42-SS	SST21-42-SS
SST13-48-SS	SST16-48-SS	SST21-48-SS

Countertop Wall Brackets

12-gauge stainless steel wall brackets for countertops.

WALL BRACKETS



Item Number	Depth	Height
CTB-30	29.750"	18"
CTB-25	24.75"	18"
CTB-24	23.75"	18"
CTB-18	17.750"	18"

Wall Mounted Sink Enclosure

Stainless steel sink enclosures are supplied with 12-gauge end support brackets, an 18-gauge apron, and 18-gauge removeable front panel.

WALL MOUNTED SINK ENCLOSURE




Item Number	Width	Depth	Height
SSE-30X23X22	30"	23"	22"
SSE-30X24X22	30"	24"	22"
SSE-30X29X22	30"	29"	22"
SSE-36X23X22	36"	23"	22"
SSE-36X24X22	36"	24"	22"
SSE-36X29X22	36"	29"	22"
SSE-48X23X22	48"	23"	22"
SSE-48X24X22	48"	24"	22"
SSE-48X29X22	48"	29"	22"

Utility Tables

All Stainless Steel Utility Tables are made with 2" square tube legs and adjustable bullet feet. Front to back and center stretchers with fully welded 4" tall aprons are standard. Adjustable legs are available and fixed heights are 35" tall, unless otherwise noted. Drawer options and lockable casters are available upon request. Adjustable height tables are manually adjustable with supplied locking clevis pins.*

UTILITY TABLES



Fixed Height



Adjustable Height


	A	B
Item Number	Length	Depth
STF-3024-36	28"	22"
STF-3030-36	28"	28"
STF-3036-36	28"	34"
STF-3624-36	34"	22"
STF-3630-36	34"	28"
STF-3636-36	34"	34"
STF-4824-36	46"	22"
STF-4830-36	46"	28"
STF-4836-36	46"	34"
STF-5424-36	52"	22"
STF-5430-36	52"	28"
STF-5436-36	52"	34"
STF-6024-36	58"	22"
STF-6030-36	58"	28"
STF-6036-36	58"	34"
STF-6624-36	64"	22"
STF-6630-36	64"	28"
STF-6636-36	64"	34"
STF-7224-36	70"	22"
STF-7230-36	70"	28"
STF-7236-36	70"	34"
STF-8424-36	82"	22"
STF-8430-36	82"	28"
STF-8436-36	82"	34"
STF-9624-36	94"	22"
STF-9630-36	94"	28"
STF-9636-36	94"	34"

*ADD "CAST" to part number for four locking casters

Table Tops

Stainless steel counter tops are 1” thick 16-gauge with stainless steel hatchanel substrate and sprayed with a sound deadening compound. Table tops are available in 24”, 30”, and 36” depths as standard. Customizations are available.

STAINLESS STEEL TABLE TOPS



	A	B	
Item Number	Length	Depth	Width
SCT-3024-1	30"	24"	1"
SCT-3624-1	36"	24"	1"
SCT-4224-1	42"	24"	1"
SCT-4824-1	48"	24"	1"
SCT-5424-1	54"	24"	1"
SCT-6024-1	60"	24"	1"
SCT-6624-1	66"	24"	1"
SCT-7224-1	72"	24"	1"
SCT-8424-1	84"	24"	1"
SCT-9624-1	96"	24"	1"



	A	B	
Item Number	Width	Depth	Height
SCT-3030-1	30"	30"	1"
SCT-3630-1	36"	30"	1"
SCT-4230-1	42"	30"	1"
SCT-4830-1	48"	30"	1"
SCT-5430-1	54"	30"	1"
SCT-6030-1	60"	30"	1"
SCT-6630-1	66"	30"	1"
SCT-7230-1	72"	30"	1"
SCT-8430-1	84"	30"	1"
SCT-9630-1	96"	30"	1"

	A	B	
Item Number	Width	Depth	Height
SCT-3036-1	30"	36"	1"
SCT-3636-1	36"	36"	1"
SCT-4236-1	42"	36"	1"
SCT-4836-1	48"	36"	1"
SCT-5436-1	54"	36"	1"
SCT-6036-1	60"	36"	1"
SCT-6636-1	66"	36"	1"
SCT-7236-1	72"	36"	1"
SCT-8436-1	84"	36"	1"
SCT-9636-1	96"	36"	1"

Scullery Sinks

Stainless steel 14-gauge scullery sinks come with 2” square tube legs on ADJ bullet feet with front to back and center stretchers. Tops are pitched for drainage. The “RFD” style comes with a lower shelf and boxed curb. Umbilical and column notches are available. Any size or configuration available. Comes with wall clips, drain outlets, and basket strainers.

SCULLERY SINKS



	A	B	
Item Number	Length	Depth	# of Sinks
SCUL-72-1B	72"	30"	1
SCUL-72-2B	72"	30"	2
SCUL-96-1B	96"	30"	2
SCUL-96-2B	96"	30"	3



Item Number	Width	Depth	Height
RFDS-72-1B	72"	33"	1
RFDS-72-2B	72"	33"	2
RFDS-96-2B	96"	33"	2
RFDS-96-3B	96"	33"	3
RFDS-120-2B	120"	33"	2
RFDS-120-3B	120"	33"	3

Canopy Hoods

All stainless steel 16-gauge canopy hoods wall mount or island style are supplied with hanging rods and attached brackets with duct collar. Standard and custom sizes available.

ISLAND CANOPY HOODS



Item Number	Width	Depth	Height
ICH-3630-30	36"	30"	30"
ICH-4830-30	48"	30"	30"
ICH-6030-30	60"	30"	30"
ICH-7230-30	72"	30"	30"
ICH-9630-30	96"	30"	30"

WALL CANOPY HOODS




Item Number	Width	Depth	Height
WCH-3630-30	36"	30"	30"
WCH-4830-30	48"	30"	30"
WCH-6030-30	60"	30"	30"
WCH-7230-30	72"	30"	30"
WCH-9630-30	96"	30"	30"

Umbilicals

Stainless steel utility chases “Umbilicals” are 18 gauge panels with 16 gauge fittings. Supplied with removeable panels and mounting collars.

ISLAND UMBILICAL



CONTACT CSI FOR PRICING AND PART NUMBERS

WALL UMBILICAL



CONTACT CSI FOR PRICING AND PART NUMBERS

Countertops

Stainless steel countertops and tables are constructed from type 304 stainless steel, polished to a #4 finish. Type 316 stainless steel is available upon request for an additional charge. Countertops and tables are constructed from 16-gauge thick material and sink bowls from 14-gauge material.

Table frames are made using 2" square tubular construction with adjustable bullet-type feet.

COUNTERTOPS



Marine Edge



Mechanical Joint

- 14-Gauge Sink Bowls
- Drain Outlets with Basket Strainers
- (OPTIONAL) Wall Mount Brackets
- Wall Clips provided
- Pitches Available
- Partial & Full End Splashes
- Tops crated for shipping
- Price includes a submittal dwg.
- Stainless steel hat channel substrate
- Sound deadening compound applied




Stainless Steel Custom Countertops


Free Standing Shelves

DISCOVERY SHELVING SYSTEM

Stainless steel shelving from type 304, polished to a #4 finish. 12-gauge parts with 14-gauge brackets available in standard and custom sizes. Table top mount and floor mount available.




Bracket Model #
DB-12-L
DB-12-R
DB-16-L
DG-16-R




Discovery Post
DP-48 2" x 3" x 48"
DP-84 2" x 3" x 84"

Wall Shelving


Available in heavy duty 12-gauge add "HD" to prefix of part number.




SB-12
Single Slot Bracket



SBD-12
Double Slot Bracket



WSB-12
Wall Mounting Bracket




SBB-12
Bookend Bracket

Single, Double, Wall Mount and Bookend shelf brackets offer support solutions for any application. All four come in 12" and 16" depths. Custom sizes are available.



Your customized workstation should have all the components to keep your materials right where you need them. That's why Canadian Scientific Lab Systems offers a wide variety of high-quality stainless steel shelving options in an array of sizes and custom options to meet your specific needs.

SS-Wall Standard SINGLE DOUBLE



These 14-gauge stainless steel standards are manufactured from one piece of steel. The sides are formed, then the top and bottom are folded to enclose the standard, eliminated an area that can collect materials from the environment.


Counter sunk screw holes leave the surface of the standard flush, so mounting hardware does not interfere with brackets or shelving.

Wall Shelving

Your customized workstation should have all the components to keep your materials right where you need them. That's why Inter Dyne Systems offers a wide variety of high-quality stainless steel shelving options in an array of sizes and custom options to meet your specific needs.




The standard wall-mounted shelf comes in a variety of sizes and can be easily mounted and moved to accommodate different materials.



All of Inter Dynes shelves are manufactured with strenght and durability in mind. The front and back edge is first bent to form the front face of the shelf, the remaining 1" of steel is bent again to give a smooth finished edge to each shelf. All shelves come in standard 12" and 16" depths. Custom sizes available.




A wall shelf with seismic lip keeps lab components from slipping off the edge of the shelf.



Seismic lip is manufactured from solid 1/4" #304 stainless steel rod with a # 4 finish. Available in two standard sizes, 11" and 14" lengths. Custom sizes available.



The shelf wrap-around option provides additional support to keep materials in place.



The shelf wrap-around option comes as a 4" standard height. Custom heights available.

Seismic Lips

All turret posts are drilled and tapped on their bottom side to accept a 1/4-20 threaded screw. Between the turrets we cut 1/4" rods to specific lengths to assure a simple and sturdy installation. Turret posts are stocked in either a middle post or a corner post with two access holes or an end post with one access hole. Canadian Scientific Lab Systems stocks these turrets and rods in order to meet any customer's time constraints.

TURRET POSTS

Model Number	Dimensions
MT-1	.75" Dia. x 1.5" H - Middle
CT-1	.75" Dia. x 1.5" H - Corner
ET-1	.75" Dia. x 1.5" H - End

Drilled and tapped 1/4-20 thread on bottom side for mounting to your shelf. 1/4" rod cut to length. 3/8" Dia. available upon request.



SEISMIC LIPS - STANDARD, STOCKED SIZES



Model Number	Dimensions
SL-11	11" L x 2" H x 1/4" Dia.
SL-14	14" L x 2" H x 1/4" Dia.
SL-15	15" L x 2" H x 1/4" Dia.
SL-11-DT	11" L x 2" H x 1/4" Dia. Drilled and tapped
SL-14-DT	14" L x 2" H x 1/4" Dia. Drilled and tapped
SL-15-DT	15" L x 2" H x 1/4" Dia. Drilled and tapped

Custom lengths available.

Acrylic

Custom sizes available, 3/4" thick is standard, 1" available upon request. All pegboards are supplied with the white pegs 6.5" long unless specified, eased and flamed polish edges. Drip troughs, bottom chamfer and u-channels available. Countersunk installation holes, as well as acrylic offset spacers available upon request.

ACRYLIC PEGBOARDS


Model Number	Dimensions	# of Pegs
APB-1824	18" x 24"	16
APB-2424	24" x 24"	25
APB-2430	24" x 30"	30
APB-243	24" x 36"	39
APB-3024	30" x 24"	33
APB-3030	30" x 30"	41
APB-3036	30" x 36"	50
APB-3624	36" x 24"	39
APB-3630	36" x 30"	50
APB-3636	36" x 36"	61
APB-4830	48" x 30"	68
APB-4836	48" X 36"	83
APB-6030	60" x 30"	33
APB-6033	60" x 33"	39
APB-6036	60" x 36"	44
APB-6633	66" x 33"	44



NOTCHED ACRYLIC SPLASH GUARD / DIVIDER PANELS

Custom sizes available. 1/2" or 3/4" thick, clear with eased and flame polished edges. 3/8" and 1/4" thickness available.


Model Number	Dimensions
NASG-2512-0.5	25" x 12" x .5"
NASG-2512-0.75	25" x 12" x .75"
NASG-3012-0.5	30" x 12" x .5"
NASG-3012-0.75	30" x 12" x .75"
NASG-2518-0.5	25" x 18" x .5"
NASG-2518-0.75	25" x 18" x .75"
NASG-3018-0.5	30" x 18" x .5"
NASG-3018-0.75	30" x 18" x .75"



ACRYLIC SPLASH GUARD / DIVIDER PANELS


Custom sizes available. 1/2" or 3/4" thick, clear with eased and flame polished edges. 3/8" and 1/4" thickness available.

Model Number	Dimensions
ASG-2512-0.5	25" x 12" x .5"
ASG-2512-0.75	25" x 12" x .75"
ASG-3012-0.5	30" x 12" x .5"
ASG-3012-0.75	30" x 12" x .75"
ASG-2518-0.5	25" x 18" x .5"
ASG-2518-0.75	25" x 18" x .75"
ASG-3018-0.5	30" x 18" x .5"
ASG-3018-0.75	30" x 18" x .75"



Drip Trough and Accessories

Stainless steel drip troughs can be mounted to the bottom of any acrylic or epoxy pegboard. All troughs come with three feet of clear PVC tubing and pre-drilled holes for easy installation. You can select a stainless steel screen insert or drain grid to create an additional shelf-like surface for drying. All sizes are stocked and ready for immediate shipping.




Stainless steel face-mount drip troughs are designed for use with any pegboard. Troughs come complete with three feet of clear PVC tubing and pre-drilled holes for easy installation. Face-mount drip troughs can be ordered with a screen insert or drain grid to add additional drying or storage space. They are manufactured from high-quality stainless steel and are stocked for immediate shipment.

PEGBOARD MOUNT




End View



FACE MOUNT



End View



Peg Board Mount	Face Mount	Screen Insert	Drain Grid
HD-18-2	FDT-18-2	SI-18-2	DG-18-2
HD-18-4	FDT-18-4	SI-18-4	DG-18-4
HD-20-2	FDT-20-2	SI-20-2	DG-20-2
HD-20-4	FDT-20-4	SI-20-4	DG-20-4
HD-24-2	FDT-24-2	SI-24-2	DG-24-2
HD-24-4	FDT-24-4	SI-24-4	DG-24-4
HD-30-2	FDT-30-2	SI-30-2	DG-30-2
HD-30-4	FDT-30-4	SI-30-4	DG-30-4
HD-32-2	FDT-32-2	SI-32-2	DG-32-2
HD-32-4	FDT-32-4	SI-32-4	DG-32-4
HD-36-2	FDT-36-2	SI-36-2	DG-36-2
HD-36-4	FDT-36-4	SI-36-4	DG-36-4
HD-42-2	FDT-42-2	SI-42-2	DG-42-2
HD-42-4	FDT-42-4	SI-42-4	DG-42-4
HD-48-2	FDT-48-2	SI-48-2	DG-48-2
HD-48-4	FDT-48-4	SI-48-4	DG-48-4

STAINLESS STEEL SCREEN INSERT



STAINLESS STEEL DRAIN GRID




The stainless steel insert & drain grid are ideal for keeping items off the bottom of the trough. They can also be utilized as a drying shelf. These items can be used with any Inter Dyne Systems drip trough or Victoria Style pegboard. 4" wide screen insert and drain grid shown at left.

Stainless Steel Pegboards

The Victoria or “V” pegboard is constructed of lightweight stainless steel, offered in a range of sizes and variety of mounting options. Along the bottom is a drip trough that funnels water away to a drain, allowing equipment to dry quickly and remain sanitized. Each “V” style pegboard includes a drip trough, wall bracket, standard 6” white pegs and 3 ft. of clear plastic drain hose.

The “B” Baron is also made of lightweight stainless steel. It is offered in three sizes and a variety of mounting options. A drip deflector at the bottom diverts water away from walls and directs it to a sink or drain. Each “B” style pegboard includes a wall bracket and standard 6” white pegs.



“V” VICTORIA STYLE “B” BARON STYLE

“V” VICTORIA STYLE PEGBOARDS



Model Number	Size W x H	# Pegs	Ship Weight
V1824	18" x 24"	15	13 lbs
V2418	24" x 18"	16	13 lbs
V2424	24" x 24"	20	16 lbs
V2430	24" x 30"	32	23 lbs
V2436	24" x 36"	40	25 lbs
V3024	30" x 24"	25	23 lbs
V3030	30" x 30"	50	24 lbs
V3036	30" x 36"	40	23 lbs
V3624	36" x 24"	30	25 lbs
V3630	36" x 30"	60	23 lbs
V3636	36" x 36"	66	25 lbs
V4824	48" x 24"	40	24 lbs
V4830	48" x 30"	48	28 lbs
V4836	48" x 36"	88	31 lbs

