

# ENTERPRISE

## Free-Standing Workstations



KEWAUNEE®

# INSPIRATION



1. Double-sided 6" square post, Snow White with Slate Kemresin®, task lighting, wood suspended cabinets.
2. Double-sided 6" round post, Snow White with Black Kemresin®, stainless shelf retaining rods, task lighting, 3 double fixtures per post, steel mobile cabinets with pullout writing top shelf and topped with Black Kemresin®.
3. Round Leg stainless steel Enterprise work stations. with stainless casework and worktops.



# INSPIRATION

- 4. 6" post single-sided with Enterprise tables. Snow White with Grey Kemresin. Heavy duty leveling casters and suspended wood on steel cabinets.
- 5. Double-sided 6" post, Snow White, with Grey Kemresin®, nickel double fixtures, ADA height workstation
- 6. Double-sided 6" round post, Snow White, with wood shelves, Slate Kemresin®, and mobile wood-on-steel casework.
- 7. Double-sided 6" post stainless steel enterprise tables with full-overlay, mobile, stainless casework.



# Free-Standing Workstations

- Self-supporting workstations
- Single- and double-sided
- Adjustable height worksurfaces
- Integrated shelving systems
- Pre-plumbed integrated services
- Pre-wired integrated power and data
- SEFA 10, Class 7, Category 4
- UL 61010-1 listed when pre-wired

A North Carolina based company founded in 1906, Kewaunee's history of quality products and scientific expertise positions it as a recognized global leader in the design, manufacture, and installation of laboratory, healthcare, and technical furniture products. Products include steel and wood casework, fume hoods, adaptable modular systems, movable workstations, stand-alone benches, and epoxy resin worksurfaces and sinks. Kewaunee's corporate headquarters are located in Statesville, North Carolina.





## Contents

|                                         |    |
|-----------------------------------------|----|
| How to Order.....                       | 6  |
| 2" Round Post Support Frames.....       | 8  |
| 6" Round Post Support Frames.....       | 12 |
| Enterprise Part Number Explanation..... | 16 |
| Electrical and Data Information.....    | 17 |
| Configured Support Frame Accs.....      | 18 |
| Blank Support Frame Accessories.....    | 20 |
| Worksurface Frames.....                 | 24 |
| Enterprise Shelving.....                | 27 |
| Free-Standing Table Frames.....         | 34 |
| Steel and Wood Cabinets.....            | 36 |
| Ceiling Service Panels.....             | 38 |
| Configured Technical Data Sheet.....    | 42 |
| Blank Frames Technical Data.....        | 48 |
| Enterprise Dimensional Data.....        | 60 |
| Alternate Chase Frame Dim Data.....     | 64 |

Kewaunee's Enterprise is a self-supporting moveable workstation table designed to support over 1,200 lb meeting SEFA 10, Adaptable Casework Systems, Class 7 standards Category 4. It is available with an above counter support structure, or as an adjustable height four legged table. Enterprise is available pre-plumbed and pre-wired and easily incorporates with overhead carrier and panel systems. Worksurfaces are adjustable from sitting height to standing height, and overhead shelving is adjustable on 1" increments. Enterprise is designed for mobile, suspended, or free-standing casework. When pre-wired, the workstation is listed to UL 61010-1.

Kewaunee's Enterprise provides a completely moveable laboratory casework system to meet the needs of a changing lab environment.

# How to Order

Enterprise Workstations are ordered through Kewaunee's **Authorized Dealer Network** who are Kewaunee's partners and have been thoroughly trained on our products.

To find your representative please visit [Kewaunee.com](http://Kewaunee.com).

The steps below will help you understand what items are needed to make up a complete workstation.

# 1

[click here to jump to this section](#)

## Select 2" or 6" Post, Round or Square Style

Select your post depth. 2" posts are for single-sided workstations only, while 6" posts are configured as double-sided.

[click here to jump to this section](#)

# 2

## Select Electrical, Data, and Services Configurations

A variety of electrical and data options are available for each bench as well as appropriate plumbing and connections for your situation.

# 3

[click here to jump to this section](#)

## Select Worksurface Frames and Worksurfaces

If your workstation is single-sided, select a single worksurface frame. Double-sided benches simply require two frames, one attached to each side of the post. Order worktops in your preferred color. One per worksurface frame.

[click here to jump to this section](#)

# 4

## Order Shelves and Any Other Accessories

Shelves and other accessories should also be selected.

# 5



[click here to jump to this section](#)

## Fill out Technical Data Sheet

Techsheets are required for each "Configuration of Workstation".


These are used to tell Kewaunee exactly where each fixture and service unit should be placed.

# Enterprise Workstation Information



**SEFA 10 Class 7**  
(on glides)

**SEFA 10 Class 8**  
(on casters)



## Single-sided Workstation

### 2" Post

[Page 8](#)

Shelves: on front only  
 Worksurface: on front only  
 Services: in post - plumbed with hard tubing  
 Electric/Data: in horizontal chase (one side only) and below worksurface with add-on box  
 Weight Rating: 1,600lb max, 200lb per shelf

## Double-sided Workstations

### 6" Post

[Page 13](#)

Shelves: on front and back  
 Worksurface: on front and back  
 Services: in post - plumbed with flexible tubing  
 Electric/Data: in horizontal chase (both sides) and below worksurface in post  
 Weight Rating: 1,600lb max per side, 200lb per shelf

# 2" Round Post Support Frames

## Specifications:


**Single-sided 2" Workstation Support Frames** consist of two, 2" diameter, powder coated steel, tubular posts with leveling feet. An upper channel and lower horizontal service chase, located 53 1/2" above the floor, are securely welded to the posts to form a rigid support frame.

The upper portion of the frame is slotted on 1" centers to accept shelving. The lower section is bored and tapped for the attachment of the adjustable height worksurface frame.

42" and 48" support frames provided with two duplex outlets; 60" and 72" support frames provided with four duplex outlets.


Frames are supplied pre-punched for electrical/data and service fixtures.

## Single-sided 2" Workstation Support Frames




**2" Round Workstation Support Frame - Configured**

48" Long **FE30C840248C**  
60" Long **FE34C840260C**  
72" Long **FE34C840272C**




**2" Round Workstation Support Frame - Configured**

60" Long **FE30C840260C**  
72" Long **FE30C840272C**



**2" Round Workstation Support Frame - Blank**

42" Long **FE30C840242U**  
48" Long **FE30C840248U**  
60" Long **FE34C840260U**  
72" Long **FE34C840272U**



**2" Round Workstation Support Frame - Blank**

60" Long **FE30C840260U**  
72" Long **FE30C840272U**

# 2" Square Post Support Frames

## Specifications:


**Single-sided 2" Workstation Support Frames** consist of two, 2", powder coated steel, tubular posts with leveling feet. An upper channel and lower horizontal service chase, located 53 1/2" above the floor, are securely welded to the posts to form a rigid support frame.

The upper portion of the frame is slotted on 1" centers to accept shelving. The lower section is bored and tapped for the attachment of the adjustable height worksurface frame.

42" and 48" support frames provided with two duplex outlets; 60" and 72" support frames provided with four duplex outlets.


Frames are supplied pre-punched for electrical/data and service fixtures.

## Single-sided 2" Workstation Support Frames




**2" Square Workstation Support Frame - Configured**

|          |                     |
|----------|---------------------|
| 48" Long | <b>FE36C840248C</b> |
| 60" Long | <b>FE37C840260C</b> |
| 72" Long | <b>FE37C840272C</b> |




**2" Square Workstation Support Frame - Configured**

|          |                     |
|----------|---------------------|
| 60" Long | <b>FE36C840260C</b> |
| 72" Long | <b>FE36C840272C</b> |



**2" Square Workstation Support Frame - Blank**

|          |                     |
|----------|---------------------|
| 42" Long | <b>FE36C840242U</b> |
| 48" Long | <b>FE36C840248U</b> |
| 60" Long | <b>FE37C840260U</b> |
| 72" Long | <b>FE37C840272U</b> |



**2" Square Workstation Support Frame -Blank**

|          |                     |
|----------|---------------------|
| 60" Long | <b>FE36C840260U</b> |
| 72" Long | <b>FE36C840272U</b> |

# 2" Blank Frames Alternate Chases


## Specifications:

**Single-sided 2" Workstation Support Frames** consist of two, 2" round/square, powder coated steel, tubular posts with leveling feet. An upper channel and lower horizontal service chase, or two chases, securely welded to the posts to form a rigid support frame. The upper portion of the frame is slotted on 1" centers to accept shelving.

The lower section is bored and tapped for the attachment of the adjustable height worksurface frame.


Frames are supplied pre-punched for electrical/data and service fixtures.

## 2" Blank Workstation Support Frames with Alternative Chase Locations




**2" Round Workstation Support Frame - Blank**

|          |                     |
|----------|---------------------|
| 42" Long | <b>FE30C840242L</b> |
| 48" Long | <b>FE30C840248L</b> |
| 60" Long | <b>FE34C840260L</b> |
| 72" Long | <b>FE34C840272L</b> |




**2" Round Workstation Support Frame - Blank**

|          |                     |
|----------|---------------------|
| 42" Long | <b>FE30C840242D</b> |
| 48" Long | <b>FE30C840248D</b> |
| 60" Long | <b>FE34C840260D</b> |
| 72" Long | <b>FE34C840272D</b> |



**2" Square Workstation Support Frame - Blank**

|          |                     |
|----------|---------------------|
| 42" Long | <b>FE36C840242L</b> |
| 48" Long | <b>FE36C840248L</b> |
| 60" Long | <b>FE37C840260L</b> |
| 72" Long | <b>FE37C840272L</b> |



**2" Square Workstation Support Frame -Blank**

|          |                     |
|----------|---------------------|
| 42" Long | <b>FE36C840242D</b> |
| 48" Long | <b>FE36C840248D</b> |
| 60" Long | <b>FE37C840260D</b> |
| 72" Long | <b>FE37C840272D</b> |

# 2" Post Configured Frame Options


## Specifications:

**Single-sided 2" Workstation Support Frames** will accept numerous service, electrical, and data/com options. Each post will accept up to three service fittings, pre-piped for easy connection to overhead services. The horizontal chase will accept both electrical receptacles and data/com outlets, and may be pre-wired for overhead connections. Each post will also accept an optional electrical fixture below the worksurface.

Workstation Support Frame Technical Data Sheets at the end of the catalog may be used to specify services and options required.

**Assembled benches are listed to UL 61010-1 when pre-wired with cataloged fixtures.**

## 2" Workstation Support Frame Plumbing and Electrical Options



# 6" Round Post Support Frames

## Specifications:

**6" Workstation Support Frames** consist of two, 2" x 6" rectangular, powder coated steel, tubular posts with leveling feet. An upper cable chase and lower horizontal service chase, located 53 1/2" above the floor, are securely fastened to the posts to form a rigid support frame.


The upper portion of the frame is slotted on both sides on 1" centers to accept shelving. The lower section is bored

and tapped for the attachment of the adjustable height worksurface.

42" and 48" support frames provided with two duplex outlets; 60" and 72" support frames provided with four duplex outlets.


Frames on this page are supplied pre-punched for electrical/data and service fixtures.

## 6" Workstation Support Frames




**6" Round Workstation Support Frame - Configured**

|          |                     |
|----------|---------------------|
| 48" Long | <b>FE30C840648C</b> |
| 60" Long | <b>FE34C840660C</b> |
| 72" Long | <b>FE34C840672C</b> |




**6" Round Workstation Support Frame - Configured**

|          |                     |
|----------|---------------------|
| 60" Long | <b>FE30C840660C</b> |
| 72" Long | <b>FE30C840672C</b> |



**6" Round Workstation Support Frame - Blank**

|          |                     |
|----------|---------------------|
| 42" Long | <b>FE30C840642U</b> |
| 48" Long | <b>FE30C840648U</b> |
| 60" Long | <b>FE34C840660U</b> |
| 72" Long | <b>FE34C840672U</b> |



**6" Round Workstation Support Frame - Blank**

|          |                     |
|----------|---------------------|
| 60" Long | <b>FE30C840660U</b> |
| 72" Long | <b>FE30C840672U</b> |

# 6" Square Post Support Frames

## Specifications:

**6" Workstation Support Frames** consist of two, 2" x 6" rectangular, powder coated steel, tubular posts with leveling feet. An upper cable chase and lower horizontal service chase, located 53 1/2" above the floor, are securely fastened to the posts to form a rigid support frame.


The upper portion of the frame is slotted on both sides on 1" centers to accept shelving. The lower section is bored

and tapped for the attachment of the adjustable height worksurface.

42" and 48" support frames provided with two duplex outlets; 60" and 72" support frames provided with four duplex outlets.


Configured frames are supplied pre-punched for electrical/ data and service fixtures.

## 6" Workstation Support Frames




**6" Square Workstation Support Frame - Configured**

|          |                     |
|----------|---------------------|
| 48" Long | <b>FE36C840648C</b> |
| 60" Long | <b>FE37C840660C</b> |
| 72" Long | <b>FE37C840672C</b> |




**6" Square Workstation Support Frame - Configured**

|          |                     |
|----------|---------------------|
| 60" Long | <b>FE36C840660C</b> |
| 72" Long | <b>FE36C840672C</b> |



**6" Square Workstation Support Frame - Blank**

|          |                     |
|----------|---------------------|
| 42" Long | <b>FE36C840642U</b> |
| 48" Long | <b>FE36C840648U</b> |
| 60" Long | <b>FE37C840660U</b> |
| 72" Long | <b>FE37C840672U</b> |



**6" Square Workstation Support Frame - Blank**

|          |                     |
|----------|---------------------|
| 60" Long | <b>FE36C840660U</b> |
| 72" Long | <b>FE36C840672U</b> |


# 6" Blank Frames Alternate Chases

## Specifications:

**6" Workstation Support Frames** consist of two, 2" x 6", powder coated steel, tubular posts with leveling feet. An upper cable chase and lower horizontal service chase, or two, securely fastened to the posts to form a rigid support frame. The upper portion of the frame is slotted on both sides on 1" centers to accept shelving.


The lower section is bored and tapped for the attachment of the adjustable height worksurface.

## 6" Blank Workstation Support Frames with Alternative Chase Locations




**6" Round Workstation Support Frame - Blank**

|          |                     |
|----------|---------------------|
| 42" Long | <b>FE30C840642L</b> |
| 48" Long | <b>FE30C840648L</b> |
| 60" Long | <b>FE34C840660L</b> |
| 72" Long | <b>FE34C840672L</b> |




**6" Round Workstation Support Frame - Blank**

|          |                     |
|----------|---------------------|
| 42" Long | <b>FE30C840642D</b> |
| 48" Long | <b>FE30C840648D</b> |
| 60" Long | <b>FE34C840660D</b> |
| 72" Long | <b>FE34C840672D</b> |



**6" Square Workstation Support Frame - Blank**

|          |                     |
|----------|---------------------|
| 42" Long | <b>FE36C840642L</b> |
| 48" Long | <b>FE36C840648L</b> |
| 60" Long | <b>FE37C840660L</b> |
| 72" Long | <b>FE37C840672L</b> |



**6" Square Workstation Support Frame -Blank**

|          |                     |
|----------|---------------------|
| 42" Long | <b>FE36C840642D</b> |
| 48" Long | <b>FE36C840648D</b> |
| 60" Long | <b>FE37C840660D</b> |
| 72" Long | <b>FE37C840672D</b> |

# 6" Post Configured Frame Options

## Specifications:


**6" Workstation Support Frames** will accept numerous service, electrical, and data/com options. Each post will accept up to three service fittings, pre-piped for easy connection to overhead services. They may be piped with hoses up to 8 feet above the top of the frame. The Horizontal Chase will accept both electrical receptacles and data/com outlets, and may be pre-wired for overhead connections.

Each post will also accept an optional electrical fixture below the worksurface.

Workstation Support Frame Technical Data Sheets, at the end of the catalog, may be used to specify services and options required.

**Assembled benches are listed to UL 61010-1 when pre-wired with cataloged fixtures.**

## 6" Blank Workstation Support Frames Plumbing and Electric



# Enterprise Part Number Explanation

## Frame Number Explanation

FE30C840672C-SKC2D3L4

### Enterprise Type

FE30 = Round Tube Enterprise Bench  
 FE34 = Round Tube No Center Divider  
 FE36 = Square Tube Enterprise Bench  
 FE37 = Square Tube No Center Divider

*FE34 & FE37 are available in 60" and 72" Length Only*

### Construction & Design

### Post Height

84 = 84" High

### Post Depth

02 = 02" Deep

06 = 06" Deep

### Required Selections

See Below for Configured Benches

### Frame Type

C = Configured

U = Upper Chase

L = Lower Chase

D = Dual Chase

*Blank  
Frames*

### Overall Length

42 = 42" Long

48 = 48" Long

60 = 60" Long

72 = 72" Long

*42 only available on blank frames*

## Configured Electrical & Data Fixtures Explanation

FE30C840672C-SKC2D3L4

### Receptacle Type

SK = Spec. Grade - Black

SV = Spec. Grade - Ivory

SW = Spec. Grade - White

SG = Spec. Grade - Grey

SR = Spec. Grade - Red

HK = Hospital Grade - Black

HV = Hospital Grade - Ivory

HW = Hospital Grade - White

HG = Hospital Grade - Grey

HR = Hospital Grade - Red

### Circuit Type *order extended cord and plug on page 10*

C1 = Single Circuit - exits Left Post

C2 = Single Circuit - exits Right Post

C3 = Double Circuit - exits Left Post

C4 = Double Circuit - exits Right Post

C5 = Double Circuit - exits Left Post

C6 = Double Circuit - exits Right Post

C7 = Multi Circuit - exits Left Post

C8 = Multi Circuit - exits Right Post

### Data Type

D0 = Cover Plates Covering Cutouts

D1 = (2) Cat6 Installed Left, Cover Plate Right

D2 = (2) Cat6 Installed Right, Cover Plate Left

D3 = (2) Cat6 Pre-Wired Left, Cover Plate Right

D4 = (2) Cat6 Pre-Wired Right, Cover Plate Left

D5 = (4) Cat6 Installed Left, Cover Plate Right

D6 = (4) Cat6 Installed Right, Cover Plate Left

D7 = (4) Cat6 Pre-Wired Left, Cover Plate Right

D8 = (4) Cat6 Pre-Wired Right, Cover Plate Left

### Leg Receptacle

L0 = No Leg Receptacle

L1 = Left Post Duplex, Box-Mount,  
Cover Plate Right

L2 = Right Post Duplex, Box-Mount,  
Cover Plate Left

L3 = Both Post Duplex, Box-Mount

L4 = Left Post Duplex, Flush-Mount

L5 = Right Post Duplex, Flush-Mount

L6 = Both Post Duplex, Flush-Mount

*Single  
Sided*


*Double  
Sided*

# Electrical and Data Information

2


## Circuit Type: C1 / C2

Requires one three-pronged power cord.




## Circuit Type: C3 / C4

Requires two three-pronged or one four-pronged power cords.




## Circuit Type: C5 / C6

Requires two three-pronged or one four-pronged power cords.



## Circuit Type: C7 / C8

Requires two four-pronged power cords.



- Circuit 1
- Circuit 2
- Circuit 3
- Circuit 4


Leg Receptacles only if specified with L1 - L6 selection.

## Configured Support Frame Accs.

### Specifications:

**Power Cords** are required to carry power from the overhead service carriers to the support frame. The SJ 12/3 and SJ 12/4 cords are available in multiple lengths. Provided plugs are Twist-Lock.

### Power Cords



#### Power Cord Part Number

**FEPCC12-03B**

#### Circuit Type

**C12** = circuit selection C1-C2 - 3-Prong Plug  
**C38** = circuit selection C3-C8 - 4-Prong Plug

#### Color

**B** = Black  
**W** = White

#### Cord Length

(above top of post)

**03** = 3'  
**04** = 4'  
**05** = 5'  
**06** = 6'  
**08** = 8'

#### Single-circuit 3-wire

120 VAC — 20A — NEMA L5-20



#### Dual-circuit 4-wire

120 VAC — 20A — NEMA L14-20



# Configured Support Frame Fixtures

2

## Specifications:

**2" Posts:** tubes are factory installed and consist of copper/ stainless steel 1/4" tubing with a female 3/8" IPS outlet escutcheon mounted in the post and a female 3/8" IPS inlet at the top of the post. If hose extension selected, the 3/8" hoses are furnished with a 3/8" IPS female outlet and a 3/8" IPS male quick-connect inlet. They are available in multiple lengths from 3 feet to 8 feet above the top of the post.

**6" Posts:** tubes are reinforced PVC/braided stainless steel. The 3/8" hoses are furnished with a 3/8" IPS female outlet and a 3/8" IPS male quick-connect inlet. They are available in multiple lengths from 3 feet to 8 feet above the top of the post.

**Internal Service Tubes** are required to mount services on the workstation support frame.

**Ball Valves** provide on/off control of gas, air, or vacuum at pressures up to 75 PSI.

**Needle Point Cocks** provide metering control of laboratory gases at pressures up to 75 PSI.

**Hose Length** is the dimension of the hose that extends above the post. The 2" post is the only workstation that allows for 0' at the top of the post, which means that there is no extension hose added. 0' extension will not receive quick connects only a threaded connection.

These fixtures are only for Configured Frames.

## Service Fittings Part Number

FE242-A1P03

### Valve Model

242 = Ball Valve 2" Post  
243 = Needle Valve 2" Post  
604 = Ball Valve 6" Post  
605 = Needle Valve 6" Post

### Service

A = Air  
G = Gas  
V = Vacuum  
N = Nitrogen  
R = Argon

### Hole Location

1 = Top  
2 = Middle  
3 = Bottom

### Hose Length

00 = 0' Ext.  
03 = 3' Ext.  
04 = 4' Ext.  
05 = 5' Ext.  
06 = 6' Ext.  
08 = 8' Ext.


### Tube/Hose Type

P = Copper or PVC  
S = Stainless


P option supplies Copper Tube for 2" post or PVC Hose for 6" post.

S option supplies Stainless Tube for 2" post or Braided Stainless Hose for 6" post.


242




243



604





605



| Service  | Fitting Color | Index Symbol | Index Color | Letter Color |
|----------|---------------|--------------|-------------|--------------|
| Air      | Orange        | AIR          | Orange      | Black        |
| Gas      | Blue          | GAS          | Blue        | White        |
| Vacuum   | Yellow        | VAC          | Yellow      | Black        |
| Nitrogen | Brown         | NIT          | Brown       | White        |
| Argon    | Violet        | AR           | Violet      | White        |

Service Fitting Colors & Index Symbols



## 2" Blank Support Frame Accessories

### Specifications:


**Service Hoses** are required to carry services from overhead service carriers to the support frame. They are available in reinforced PVC or braided stainless steel.

**Power Cords** are required to carry power from the overhead service carriers to the support frame. The SJ 12/3 and SJ 12/4 cords are available in multiple lengths.

**Keyed Quick-Connect Fittings** are available for hoses


and internal service tubes to ease service connection and disconnection during bench relocation. They are available in eight sizes (keys) which keeps service outlets from being connected to the wrong service.

### Extension Service Hoses



Specify Hose and Power Cord locations on *Workstation Support Frame Technical Data Sheet*

### Power Cords



W-5621, W-5622, and W-5624 parts require a color selection of Black (**BLK**) or White (**WHT**) when ordered. **Ex: W-5621-36-BLK**

### Reinforced PVC Hose

for Non-burning Gases

3/8" IPS Female Inlet — 3/8" IPS Male Outlet

| Length | Part Number  | Color      |
|--------|--------------|------------|
| 36"    | W-1133-36-H_ | 2 = Orange |
| 48"    | W-1133-48-H_ | 4 = Yellow |
| 60"    | W-1133-60-H_ | 5 = Blue   |
| 72"    | W-1133-72-H_ | 6 = White  |
| 96"    | W-1133-96-H_ | 7 = Purple |
|        |              | 8 = Brown  |

### Braided Stainless Steel Hose

for Natural & Reactive Gases

3/8" IPS Female Inlet — 3/8" IPS Male Outlet

| Length | Part Number  |
|--------|--------------|
| 36"    | W-1233-36-H0 |
| 48"    | W-1233-48-H0 |
| 60"    | W-1233-60-H0 |
| 72"    | W-1233-72-H0 |
| 96"    | W-1233-96-H0 |

### Single-circuit 3-wire

120 VAC — 20A — NEMA 5-20

| Part Number | Length  | Plug Type      |
|-------------|---------|----------------|
| W-5621-36-  | 36"ATOP | Straight Blade |
| W-5621-48-  | 48"ATOP | Straight Blade |
| W-5621-60-  | 60"ATOP | Straight Blade |
| W-5621-72-  | 72"ATOP | Straight Blade |
| W-5621-96-  | 96"ATOP | Straight Blade |

120 VAC — 20A — NEMA L5-20

| Part Number | Length  | Plug Type  |
|-------------|---------|------------|
| W-5622-36-  | 36"ATOP | Twist-Lock |
| W-5622-48-  | 48"ATOP | Twist-Lock |
| W-5622-60-  | 60"ATOP | Twist-Lock |
| W-5622-72-  | 72"ATOP | Twist-Lock |
| W-5622-96-  | 96"ATOP | Twist-Lock |

### Dual-circuit 4-wire

120 VAC — 20A — NEMA L14-20


| Part Number | Length  | Plug Type  |
|-------------|---------|------------|
| W-5624-36-  | 36"ATOP | Twist-Lock |
| W-5624-48-  | 48"ATOP | Twist-Lock |
| W-5624-60-  | 60"ATOP | Twist-Lock |
| W-5624-72-  | 72"ATOP | Twist-Lock |
| W-5624-96-  | 96"ATOP | Twist-Lock |

### Keyed Quick-Connect Fittings for Hoses



**Male  
Hose Mount**  
3/8" IPS Male

| Part Number | Key Color |
|-------------|-----------|
| W-8111-K2   | Orange    |
| W-8111-K4   | Yellow    |
| W-8111-K5   | Blue      |
| W-8111-K7   | Purple    |
| W-8111-K8   | Brown     |



**Female  
Hose Mount**  
3/8" IPS Female

| Part Number | Key Color |
|-------------|-----------|
| W-9112-K2   | Orange    |
| W-9112-K4   | Yellow    |
| W-9112-K5   | Blue      |
| W-9112-K7   | Purple    |
| W-9112-K8   | Brown     |

| Key Number and Color Code | Standard Service | Female Body Sleeve Dia. |      | Male Body Sleeve Dia. |      |
|---------------------------|------------------|-------------------------|------|-----------------------|------|
|                           |                  | in.                     | mm   | in.                   | mm   |
| K2 Orange                 | Air              | 1.16                    | 29.5 | 1.02                  | 25.8 |
| K4 Yellow                 | Vacuum           | 1.22                    | 31.0 | 1.08                  | 27.4 |
| K5 Blue                   | Gas              | 1.24                    | 31.5 | 1.11                  | 28.1 |
| K7 Purple                 | Argon            | 1.31                    | 33.3 | 1.17                  | 29.7 |
| K8 Brown                  | Nitrogen         | 1.34                    | 34.0 | 1.20                  | 30.4 |

# 2" Blank Support Frame Accessories

2

## Specifications:

**Internal Service Tubes** are required to mount services on blank support frames. They are factory installed and consist of either copper or stainless steel 1/4" tubing with a female 3/8" IPS outlet escutcheon mounted in the 2" post and a female 3/8" IPS inlet at the top of the post. They are available in three lengths that correspond to the three mounting locations in the post.

**Ball Valves** are chrome plated brass and provide on/off control of gas, air, or vacuum at pressures up to 75 PSI.

**Needle Point Cocks** are chrome plated brass with stainless steel floating cone and removable brass seat. They provide metering control of laboratory gases at pressures up to 75 PSI.

## Internal Service Tubes



### Copper Tube

for Non-burning Gases

3/8" IPS Female Inlet — 3/8" IPS Female Outlet

LengthPart Number

|     |                  |             |        |
|-----|------------------|-------------|--------|
| 40" | <b>W-2321-40</b> | Top Hole    | Round  |
| 44" | <b>W-2321-44</b> | Middle Hole |        |
| 48" | <b>W-2321-48</b> | Bottom Hole |        |
| 40" | <b>W-3320-40</b> | Top Hole    | Square |
| 44" | <b>W-3320-44</b> | Middle Hole |        |
| 48" | <b>W-3320-48</b> | Bottom Hole |        |

### Stainless Steel Tube


for Natural & Reactive Gases

3/8" IPS Female Inlet — 3/8" IPS Female Outlet

LengthPart Number

|     |                  |             |        |
|-----|------------------|-------------|--------|
| 40" | <b>W-2421-40</b> | Top Hole    | Round  |
| 44" | <b>W-2421-44</b> | Middle Hole |        |
| 48" | <b>W-2421-48</b> | Bottom Hole |        |
| 40" | <b>W-3420-40</b> | Top Hole    | Square |
| 44" | <b>W-3420-44</b> | Middle Hole |        |
| 48" | <b>W-3420-48</b> | Bottom Hole |        |


## Service Fittings



### Ball Valve

3/8" IPS Male Inlet

| Part Number      | Service  |
|------------------|----------|
| <b>W-0242-0A</b> | Air      |
| <b>W-0242-0G</b> | Gas      |
| <b>W-0242-0N</b> | Nitrogen |
| <b>W-0242-0R</b> | Argon    |
| <b>W-0242-0V</b> | Vacuum   |




### Needle Point Cock

3/8" IPS Male Inlet

| Part Number      | Service  |
|------------------|----------|
| <b>W-0243-0A</b> | Air      |
| <b>W-0243-0G</b> | Gas      |
| <b>W-0243-0N</b> | Nitrogen |
| <b>W-0243-0R</b> | Argon    |
| <b>W-0243-0V</b> | Vacuum   |

## Power Bar Electrical Fixtures

Specify Power Bar Electrical Fixture locations on Workstation Support Technical Data Sheet



**0720-0K-6A**

### Type

**0720** = Data RJ45  
**S581** = Power 20 amp Spec  
**H581** = Power 20 amp Hospital  
**S656** = Power GFCI 20 amp Spec

### Color

**K** = Black  
**W** = White  
**G** = Grey  
**V** = Ivory  
**R** = Red


### Data Category (Data Only)

**5e** = Category 5e  
**6** = Category 6  
**6a** = Category 6a

## Rear Support Post Electrical Fixtures

**2" Square Frame Lower Electrical** is integrated into the post, and does not require an external junction box. Order using part number from Power Bar Electrical Section to the left.

Specify Post Electrical Fixture locations on Workstation Support Frame Technical Data Sheet



### Color

**K** = Black  
**W** = White  
**G** = Grey  
**V** = Ivory  
**R** = Red

### 120 VAC Outlet

| Part Number           | Device             |
|-----------------------|--------------------|
| <b>FE01D050202-__</b> | 20 amp Duplex      |
| <b>FE01G050202-__</b> | 20 amp GFCI Duplex |

## 6" Blank Support Frame Accessories


### Specifications:

**Internal Service Hoses** are required to mount services on the workstation support frame. They are factory installed when ordered with the service fitting and are available in either reinforced PVC or braided stainless steel. The  $\frac{3}{8}$ " hoses are furnished with a  $\frac{3}{8}$ " IPS female outlet and either a  $\frac{3}{8}$ " IPS female inlet or  $\frac{3}{8}$ " IPS male quick-connect inlet. They are available in multiple lengths from 3 feet to 8 feet above the top of the post.

**Double Ball Valves** consist of a 90° wye turret with two chrome plated valves which provide on/off control of gas, air, or vacuum at pressures up to 75 PSI.

**Double Needle Point Cocks** consist of a 90° wye turret with two chrome plated valves with stainless steel floating cone and removable brass seat. They provide metering control of laboratory gases at pressures up to 75 PSI.

### Internal Service Hoses



#### Color

- 2 = Orange
- 4 = Yellow
- 5 = Blue
- 6 = White
- 7 = Purple
- 8 = Brown

More selections necessary inside Connect to create a final part number.

Example:  
W-6136E36-H6K4V

#### Reinforced PVC Hose

for Non-burning Gases  
 $\frac{3}{8}$ " IPS Female Inlet  
 $\frac{3}{8}$ " IPS Female Outlet

| Part Number  | Length (Above Top Of Post) |
|--------------|----------------------------|
| W-6134E36-H_ | 36" ATOP                   |
| W-6134E48-H_ | 48" ATOP                   |
| W-6134E60-H_ | 60" ATOP                   |
| W-6134E72-H_ | 72" ATOP                   |
| W-6134E96-H_ | 96" ATOP                   |

#### Reinforced PVC Hose with Quick-connect Fitting

for Non-burning Gases  
 $\frac{3}{8}$ " IPS Male Quick-connect Inlet (not illustrated)  
 $\frac{3}{8}$ " IPS Female Outlet

| Part Number    | Length (Above Top Of Post) |
|----------------|----------------------------|
| W-6136E36-H_K_ | 36" ATOP                   |
| W-6136E48-H_K_ | 48" ATOP                   |
| W-6136E60-H_K_ | 60" ATOP                   |
| W-6136E72-H_K_ | 72" ATOP                   |
| W-6136E96-H_K_ | 96" ATOP                   |

#### Braided Stainless Steel Hose

for Natural & Reactive Gases  
 $\frac{3}{8}$ " IPS Female Inlet  
 $\frac{3}{8}$ " IPS Female Outlet


| Part Number  | Length (Above Top Of Post) |
|--------------|----------------------------|
| W-6234E36-HO | 36" ATOP                   |
| W-6234E48-HO | 48" ATOP                   |
| W-6234E60-HO | 60" ATOP                   |
| W-6234E72-HO | 72" ATOP                   |
| W-6234E96-HO | 96" ATOP                   |

#### Braided Stainless Steel Hose with Quick-connect Fitting

for Natural & Reactive Gases  
 $\frac{3}{8}$ " IPS Male Quick-connect Inlet (not illustrated)  
 $\frac{3}{8}$ " IPS Female Outlet

| Part Number    | Length (Above Top Of Post) |
|----------------|----------------------------|
| W-6236E36-HOK_ | 36" ATOP                   |
| W-6236E48-HOK_ | 48" ATOP                   |
| W-6236E60-HOK_ | 60" ATOP                   |
| W-6236E72-HOK_ | 72" ATOP                   |
| W-6236E96-HOK_ | 96" ATOP                   |


### Service Fittings



#### Double Ball Valve

$\frac{3}{8}$ " IPS Female Inlet

| Part Number | Service  |
|-------------|----------|
| W-0604-2A   | Air      |
| W-0604-2G   | Gas      |
| W-0604-2N   | Nitrogen |
| W-0604-2R   | Argon    |
| W-0604-2V   | Vacuum   |




#### Double Needle Point Cock

$\frac{3}{8}$ " IPS Female Inlet

| Part Number | Service  |
|-------------|----------|
| W-0605-2A   | Air      |
| W-0605-2G   | Gas      |
| W-0605-2N   | Nitrogen |
| W-0605-2R   | Argon    |
| W-0605-2V   | Vacuum   |

Specify Service Fitting locations on *Workstation Support Frame Technical Data Sheet*

| Key Number and Color Code | Standard Service | Female Body Sleeve Dia. |      | Male Body Sleeve Dia. |      |
|---------------------------|------------------|-------------------------|------|-----------------------|------|
|                           |                  | in.                     | mm   | in.                   | mm   |
| K2 Orange                 | Air              | 1.16                    | 29.5 | 1.02                  | 25.8 |
| K4 Yellow                 | Vacuum           | 1.22                    | 31.0 | 1.08                  | 27.4 |
| K5 Blue                   | Gas              | 1.24                    | 31.5 | 1.11                  | 28.1 |
| K7 Purple                 | Argon            | 1.31                    | 33.3 | 1.17                  | 29.7 |
| K8 Brown                  | Nitrogen         | 1.34                    | 34.0 | 1.20                  | 30.4 |



# 6" Blank Support Frame Accessories


2

## Specifications:

**Power Cords** are required to carry power from the overhead service carriers to the support frame. The SJ 12/3 and SJ 12/4 cords are available in multiple lengths.


### Power Bar Electrical Fixtures

Specify Power Bar Electrical Fixture locations on Workstation Support Technical Data Sheet



### Support Post Electrical Fixtures

Specify Post Electrical Fixture locations on Workstation Support Frame Technical Data Sheet



**0720-0K-6A**

#### Type

**0720** = Data RJ45  
**S581** = Power 20 amp Spec  
**H581** = Power 20 amp Hospital  
**S656** = Power GFCI 20 amp Spec


#### Color

**K** = Black  
**W** = White  
**G** = Grey  
**V** = Ivory  
**R** = Red

#### Data Category (Data Only)

**5e** = Category 5e  
**6** = Category 6  
**6a** = Category 6a

### Power Cords



W-5621, W-5622, and W-5624 parts require a color selection of Black (**BLK**) or White (**WHT**) when ordered. **Ex: W-5621-36-BLK**

#### Single-circuit 3-wire

120 VAC — 20A — NEMA 5-20

| Part Number       | Length  | Plug Type      |
|-------------------|---------|----------------|
| <b>W-5621-36-</b> | 36"ATOP | Straight Blade |
| <b>W-5621-48-</b> | 48"ATOP | Straight Blade |
| <b>W-5621-60-</b> | 60"ATOP | Straight Blade |
| <b>W-5621-72-</b> | 72"ATOP | Straight Blade |
| <b>W-5621-96-</b> | 96"ATOP | Straight Blade |

120 VAC — 20A — NEMA L5-20

| Part Number       | Length  | Plug Type  |
|-------------------|---------|------------|
| <b>W-5622-36-</b> | 36"ATOP | Twist-Lock |
| <b>W-5622-48-</b> | 48"ATOP | Twist-Lock |
| <b>W-5622-60-</b> | 60"ATOP | Twist-Lock |
| <b>W-5622-72-</b> | 72"ATOP | Twist-Lock |
| <b>W-5622-96-</b> | 96"ATOP | Twist-Lock |

#### Dual-circuit 4-wire

120 VAC — 20A — NEMA L14-20

| Part Number       | Length  | Plug Type  |
|-------------------|---------|------------|
| <b>W-5624-36-</b> | 36"ATOP | Twist-Lock |
| <b>W-5624-48-</b> | 48"ATOP | Twist-Lock |
| <b>W-5624-60-</b> | 60"ATOP | Twist-Lock |
| <b>W-5624-72-</b> | 72"ATOP | Twist-Lock |
| <b>W-5624-96-</b> | 96"ATOP | Twist-Lock |


Specify Power Cord locations on Workstation Support Frame Technical Data Sheet

# 3 Worksurface Frames

## Specifications:

**Worksurface Frames** are available in four lengths to match the workstation support frames, and are adjustable in height from 29" - 36" on 1" increments. The frames are fabricated from powder coated steel with 2" diameter or 2" square, telescoping legs. They are furnished with required installation hardware.

### Round Worksurface Frames



#### 22" Deep Worktop Frame for use with 24" Worktops

|       |              |
|-------|--------------|
| 42" L | FE80C362242N |
| 48" L | FE80C362248N |
| 60" L | FE80C362260N |
| 72" L | FE80C362272N |

#### 25" Deep Worktop Frame for use with 27" Worktops

|       |              |
|-------|--------------|
| 42" L | FE80C362542N |
| 48" L | FE80C362548N |
| 60" L | FE80C362560N |
| 72" L | FE80C362572N |

#### 28" Deep Worktop Frame for use with 30" Worktops

|       |              |
|-------|--------------|
| 42" L | FE80C362842N |
| 48" L | FE80C362848N |
| 60" L | FE80C362860N |
| 72" L | FE80C362872N |


#### 31" Deep Worktop Frame for use with 33" Worktops

|       |              |
|-------|--------------|
| 42" L | FE80C363142N |
| 48" L | FE80C363148N |
| 60" L | FE80C363160N |
| 72" L | FE80C363172N |

#### 34" Deep Worktop Frame for use with 36" Worktops

|       |              |
|-------|--------------|
| 42" L | FE80C363442N |
| 48" L | FE80C363448N |
| 60" L | FE80C363460N |
| 72" L | FE80C363472N |

### Square Worksurface Frames



#### 22" Deep Worktop Frame for use with 24" Worktops

|       |              |
|-------|--------------|
| 42" L | FE82C362242N |
| 48" L | FE82C362248N |
| 60" L | FE82C362260N |
| 72" L | FE82C362272N |

#### 25" Deep Worktop Frame for use with 27" Worktops

|       |              |
|-------|--------------|
| 42" L | FE82C362542N |
| 48" L | FE82C362548N |
| 60" L | FE82C362560N |
| 72" L | FE82C362572N |

#### 28" Deep Worktop Frame for use with 30" Worktops

|       |              |
|-------|--------------|
| 42" L | FE82C362842N |
| 48" L | FE82C362848N |
| 60" L | FE82C362860N |
| 72" L | FE82C362872N |

#### 31" Deep Worktop Frame for use with 33" Worktops

|       |              |
|-------|--------------|
| 42" L | FE82C363142N |
| 48" L | FE82C363148N |
| 60" L | FE82C363160N |
| 72" L | FE82C363172N |

#### 34" Deep Worktop Frame for use with 36" Worktops



|       |              |
|-------|--------------|
| 42" L | FE82C363442N |
| 48" L | FE82C363448N |
| 60" L | FE82C363460N |
| 72" L | FE82C363472N |

# Worksurface Frame Worktops

## Specifications:

**Worksurface Frame Worktops** are available in Kemresin Epoxy Resin or Phenolic Resin. All are 1" thick, have a 1" front overhang, a 3" or 6" rear overhang, and are flush with the ends of the worksurface frames. Rear corners are notched to fit around workstation support frames. All edges are finished.

## Worksurface Worktops with 3" Rear Notch (for use with 2" Single-sided and 6" Double-sided Support Frames)



## Worktop for Worksurface Frames Part Number

**FE95R012442N-BK**

### Top Configuration

- 95** = Round Tube
- 92** = Square Tube
- 97** = Round 6" Rear Notch
- 99** = Square 6" Rear Notch

### Material

- R** = Kemresin Epoxy
- T** = Phenolic Resin

### Worktop Thickness

- 01** = 1"

### Color

- BK** = Black
- GR** = Grey
- SL** = Slate
- PT** = Putty
- ABK** = \*Black

### Worktop Length

- 42** = 42"
- 48** = 48"
- 60** = 60"
- 72** = 72"

### Worktop Depth

- 24** = 24"
- 27** = 27"
- 30** = 30"
- 33** = 33"
- 36** = 36"


\***ABK** color is for use with Phenolic Resin only.

**Worktop Depths** are measured two inches deeper than their corresponding Worksurface Frames.

24" deep worktops for 22" frames  
27" deep worktops for 25" frames  
30" deep worktops for 28" frames  
33" deep worktops for 31" frames  
36" deep worktops for 34" frames

**FE97** and **FE99** tops provide an extended rear overhang for single-sided, 6" post benches allowing the worktop to be flush at the rear of the posts.

**Worktop to Column Dimensions** are available in the Dimensions section of this Catalog.  
[CLICK HERE TO VIEW](#)



# 3 Worksurface Frame Accessories

## Specifications:

**Rear Adjustable Legs** match the telescoping legs at the front of worksurface frames. They may be attached to any of the worksurface frames to create a free-standing table frame. (sold in pairs)


**Post Brackets** are designed to clamp two FE30 single-sided, 2" post workstations together, back-to-back, to form a double-sided assembly. (brackets are sold in pairs)

**Cabinet Hanging Rail Kit** may be added to the 28", 31", or 34" deep worksurface frame to hang suspended base cabinets.

**Cabinet Stops** are designed to locate mobile base cabinets in line with the front of the workstation. Fabricated from 16 gauge, powder coated, steel, they are furnished with all required installation hardware.

**Powder Coated Steel Back Splash** is 5" high and mounted to the worksurface frame to provide a 4" high splash. Fabricated from 11 gauge powder coated steel, it is furnished complete with installation hardware. Available in black only.

### Rear Adjustable Legs




30" - 36"H  
30" - 36"H

FE88M360202P  
FE87C360202P

round pair  
square pair


### Cabinet Hanging Rail Kit



42" L FE89M020242N  
48" L FE89M020248N  
60" L FE89M020260N  
72" L FE89M020272N

For 28", 31", 34" Deep Worksurface Frames

### Cabinet Stops



#### 25" Deep Frame Cabinet Stops

Inset Cabinets

42" L FE85C022542N  
48" L FE85C022548N  
60" L FE85C022560N  
72" L FE85C022572N

Overlay Cabinets


42" L FE85C022542V  
48" L FE85C022548V  
60" L FE85C022560V  
72" L FE85C022572V

#### 28", 31", 34" Deep Frame Cabinet Stops

Inset or Overlay

42" L FE85C020042N  
48" L FE85C020048N  
60" L FE85C020060N  
72" L FE85C020072N


### Round Tube Post Brackets



FE02M020206P pair

each bracket includes two (2) half brackets with nuts and bolts for assembly.

### Black Powder Coated Steel Back Splash




42" L FE95M040342N  
48" L FE95M040348N  
60" L FE95M040360N  
72" L FE95M040372N

# Enterprise Shelving

**Features:**

- May be used with Workstation Support Frames or individual Wall Standards
- Shelf Retaining Lips available for front or back edge of all shelves
- Task Lighting available for mounting under shelves
- Wall Standards available in heights of 30", 36", 42", 48" and 84"
- Shelving available in lengths of 29", 35", 40", 46" and 47" to match Workstation Support Frames
- Shelving available with or without front retaining rods



# 4 Wall Standards & Shelving Acc.

## Specifications:

**Enterprise Round Wall Standards** are fabricated from 11 gauge, 2" diameter powder coated steel tubing, with 1" steel stand-offs for mounting on wall. Full height wall standard is furnished with leveling foot.


All Round Enterprise Wall Standards are slotted on 1" centers and are designed to accept two shelf brackets mounted side-by-side.

**Wall mounting hardware is not included**

## Wall Standards



**Full Height Wall Standard**  
84" H    **FE10M840202N**    each



### Upper Wall Standards

|       |                     |      |
|-------|---------------------|------|
| 30" L | <b>FE10M300202N</b> | each |
| 36" L | <b>FE10M360202N</b> | each |
| 42" L | <b>FE10M420202N</b> | each |
| 48" L | <b>FE10M480202N</b> | each |

# Adjustable Height Shelves

## Specifications:


**Adjustable Height Shelves** consist of a pair of powder coated steel shelf brackets and a shelf. Shelves are adjustable in height on 1" centers from 55" to 80" above floor. Shelves overhang the rear of the shelf bracket by 3" and are available in powder coated steel, oak or maple veneer plywood.

**Powder Coated Steel Shelves** — 1" thick, formed from cold rolled steel.


**Wood Shelves** — 1" thick hardwood plywood with matching veneer on all exposed surfaces.

Shelf Lengths listed below are nominal and include 1" on either side accounting for the brackets.

## Adjustable Shelves



*Top Shelf Assembly illustrated on 2" post. Same shelf fits 6" post.*



### Frames with No Center Divider

|           |         |                   |
|-----------|---------|-------------------|
| 42" Shelf | used on | 42" Support Frame |
| 48" Shelf | used on | 48" Support Frame |
| 60" Shelf | used on | 60" Support Frame |
| 72" Shelf | used on | 72" Support Frame |

### Frames with Center Divider

|           |         |                   |
|-----------|---------|-------------------|
| 30" Shelf | used on | 60" Support Frame |
| 36" Shelf | used on | 72" Support Frame |

### 9" Deep Bracket – 12" Deep Shelf

|       | Powder Coated Steel | Oak Veneer   | Maple Veneer |
|-------|---------------------|--------------|--------------|
| 30" L | FE21C040930P        | FE21W040930R | FE21W040930M |
| 36" L | FE21C040936P        | FE21W040936R | FE21W040936M |
| 42" L | FE21C040942P        | FE21W040942R | FE21W040942M |
| 48" L | FE21C040948P        | FE21W040948R | FE21W040948M |
| 60" L | FE21C040960P        |              |              |
| 72" L | FE21C040972P        |              |              |

### 12" Deep Bracket – 15" Deep Shelf

|       | Powder Coated Steel | Oak Veneer   | Maple Veneer |
|-------|---------------------|--------------|--------------|
| 30" L | FE21C051230P        | FE21W051230R | FE21W051230M |
| 36" L | FE21C051236P        | FE21W051236R | FE21W051236M |
| 42" L | FE21C051242P        | FE21W051242R | FE21W051242M |
| 48" L | FE21C051248P        | FE21W051248R | FE21W051248M |
| 60" L | FE21C051260P        |              |              |
| 72" L | FE21C051272P        |              |              |

### 18" Deep Bracket – 21" Deep Shelf

|       | Powder Coated Steel | Oak Veneer   | Maple Veneer |
|-------|---------------------|--------------|--------------|
| 30" L | FE21C061830P        | FE21W061830R | FE21W061830M |
| 36" L | FE21C061836P        | FE21W061836R | FE21W061836M |
| 42" L | FE21C061842P        | FE21W061842R | FE21W061842M |
| 48" L | FE21C061848P        | FE21W061848R | FE21W061848M |
| 60" L | FE21C061860P        |              |              |
| 72" L | FE21C061872P        |              |              |

# 4 Adjustable Height Shelves

## Specifications:


**Adjustable Height Shelves with Retaining Rod** consist of a pair of powder coated steel shelf brackets, a shelf, and a  $\frac{3}{8}$ " diameter powder coated steel rod fastened to the front of the shelf brackets. Shelves are adjustable in height on 1" centers from 55" to 80" above floor. Shelves overhang the rear of the shelf bracket by 3" and are available in powder coated steel, oak or maple veneer plywood.

**Powder Coated Steel Shelves** — 1" thick, formed from cold rolled steel.

**Wood Shelves** — 1" thick hardwood plywood with matching veneer on all exposed surfaces.

Shelf Lengths listed below are nominal and include 1" on either side accounting for the brackets.


## Adjustable Shelves — with Retaining Rod



Frames with No Center Divider

|           |         |                   |
|-----------|---------|-------------------|
| 42" Shelf | used on | 42" Support Frame |
| 48" Shelf | used on | 48" Support Frame |
| 60" Shelf | used on | 60" Support Frame |
| 72" Shelf | used on | 72" Support Frame |

*Top Shelf Assembly illustrated on 2" post. Same shelf fits 6" post.*



Frames with Center Divider

|           |         |                   |
|-----------|---------|-------------------|
| 30" Shelf | used on | 60" Support Frame |
| 36" Shelf | used on | 72" Support Frame |

### 9" Deep Bracket - 12" Deep Shelf

|       | Powder Coated Steel | Oak Veneer   | Maple Veneer |
|-------|---------------------|--------------|--------------|
| 30" L | FE22C040930P        | FE22W040930R | FE22W040930M |
| 36" L | FE22C040936P        | FE22W040936R | FE22W040936M |
| 42" L | FE22C040942P        | FE22W040942R | FE22W040942M |
| 48" L | FE22C040948P        | FE22W040948R | FE22W040948M |
| 60" L | FE22C040960P        |              |              |
| 72" L | FE22C040972P        |              |              |

### 12" Deep Bracket - 15" Deep Shelf

|       | Powder Coated Steel | Oak Veneer   | Maple Veneer |
|-------|---------------------|--------------|--------------|
| 30" L | FE22C051230P        | FE22W051230R | FE22W051230M |
| 36" L | FE22C051236P        | FE22W051236R | FE22W051236M |
| 42" L | FE22C051242P        | FE22W051242R | FE22W051242M |
| 48" L | FE22C051248P        | FE22W051248R | FE22W051248M |
| 60" L | FE22C051260P        |              |              |
| 72" L | FE22C051272P        |              |              |

### 18" Deep Bracket - 21" Deep Shelf

|       | Powder Coated Steel | Oak Veneer   | Maple Veneer |
|-------|---------------------|--------------|--------------|
| 30" L | FE22C061830P        | FE22W061830R | FE22W061830M |
| 36" L | FE22C061836P        | FE22W061836R | FE22W061836M |
| 42" L | FE22C061842P        | FE22W061842R | FE22W061842M |
| 48" L | FE22C061848P        | FE22W061848R | FE22W061848M |
| 60" L | FE22C061860P        |              |              |
| 72" L | FE22C061872P        |              |              |

## Top Shelves

### Specifications:


**Single-sided Support Frame Top Shelves** overhang the rear of the shelf bracket by 1" and are available in powder coated steel, oak or maple veneer plywood. Single-sided support frame top shelves consist of a pair of powder coated steel shelf brackets and a shelf. Top Shelves are designed to mount at the top of the support assembly but are adjustable in height on 1" centers from 60" to 84" above floor.

**Powder Coated Steel Shelves** — 1" thick, formed from cold rolled steel.


**Wood Shelves** — 1" thick hardwood plywood with matching veneer on all exposed surfaces.

Shelf Lengths listed below are nominal and include 1" on either side accounting for the brackets.

### Top Shelf Assembly for 2" & 6" Post



*Top Shelf Assembly illustrated on 2" post. Same shelf fits 6" post.*



#### Frames with No Center Divider

|           |         |                   |
|-----------|---------|-------------------|
| 42" Shelf | used on | 42" Support Frame |
| 48" Shelf | used on | 48" Support Frame |
| 60" Shelf | used on | 60" Support Frame |
| 72" Shelf | used on | 72" Support Frame |

#### Frames with Center Divider

|           |         |                   |
|-----------|---------|-------------------|
| 30" Shelf | used on | 60" Support Frame |
| 36" Shelf | used on | 72" Support Frame |

### 9" Deep Bracket - 10" Deep Shelf

|       | Powder Coated Steel | Oak Veneer   | Maple Veneer |
|-------|---------------------|--------------|--------------|
| 30" L | FE28C040930P        | FE28W040930R | FE28W040930M |
| 36" L | FE28C040936P        | FE28W040936R | FE28W040936M |
| 42" L | FE28C040942P        | FE28W040942R | FE28W040942M |
| 48" L | FE28C040948P        | FE28W040948R | FE28W040948M |
| 60" L | FE28C040960P        |              |              |
| 72" L | FE28C040972P        |              |              |

### 12" Deep Bracket - 13" Deep Shelf

|       | Powder Coated Steel | Oak Veneer   | Maple Veneer |
|-------|---------------------|--------------|--------------|
| 30" L | FE28C051230P        | FE28W051230R | FE28W051230M |
| 36" L | FE28C051236P        | FE28W051236R | FE28W051236M |
| 42" L | FE28C051242P        | FE28W051242R | FE28W051242M |
| 48" L | FE28C051248P        | FE28W051248R | FE28W051248M |
| 60" L | FE28C051260P        |              |              |
| 72" L | FE28C051272P        |              |              |

### 18" Deep Bracket - 19" Deep Shelf

|       | Powder Coated Steel | Oak Veneer   | Maple Veneer |
|-------|---------------------|--------------|--------------|
| 30" L | FE28C061830P        | FE28W061830R | FE28W061830M |
| 36" L | FE28C061836P        | FE28W061836R | FE28W061836M |
| 42" L | FE28C061842P        | FE28W061842R | FE28W061842M |
| 48" L | FE28C061848P        | FE28W061848R | FE28W061848M |
| 60" L | FE28C061860P        |              |              |
| 72" L | FE28C061872P        |              |              |

## Square Wall Standards

### Specifications:


Enterprise Square Wall Standards are fabricated from 12 gauge powder coated steel.

All square Enterprise Wall Standards are slotted on 1" centers and are designed to accept two shelf brackets mounted side-by-side.

Wall mounting hardware is not included


Ordered each

### Square Wall Standards



#### 2" Wide Single Slot, Square Wall Standard


|       |              |
|-------|--------------|
| 30" H | FE12C300102N |
| 36" H | FE12C360102N |
| 42" H | FE12C420102N |
| 48" H | FE12C480102N |
| 84" H | FE12C840102N |



#### 4" Wide Double Slot, Square Wall Standard

|       |              |
|-------|--------------|
| 30" H | FE12C300104N |
| 36" H | FE12C360104N |
| 42" H | FE12C420104N |
| 48" H | FE12C480104N |
| 84" H | FE12C840104N |

### Shelving for Square Wall Standards



### Specifications:

Shelving for Square Wall Standards consist of a pair of powder coated steel shelf brackets and a shelf. Shelves are adjustable in height on 1" centers from 55" to 80" above floor. Shelves overhang the rear of the shelf bracket by 1".

FE23C040930P

#### Valve Model

23C = Shelf  
24C = Shelf w/ Retaining Rod

#### Shelf Depth

0409 = 10" Depth  
0512 = 13" Depth  
0618 = 19" Depth

#### Shelf Length

30 = 30"  
36 = 36"  
42 = 42"  
48 = 48"  
60 = 60"  
72 = 72"

# Shelving Accessories


## Specifications:

**Adjustable Height Shelf Brackets and Top Shelf Brackets** are for custom applications or for replacements. These are the same brackets that are used on Wood Veneer shelf assemblies. Shelf Brackets are furnished in pairs.

**Shelf Retaining Lips** are fabricated from 16 gauge powder coated steel, and may be added to the front or back of any Enterprise shelf. (required hardware included)


**Part Number is for single lip.**

### Top Shelf Brackets (in pairs)




| Depth | Part Number         |      |
|-------|---------------------|------|
| 9"    | <b>FE27C040902P</b> | pair |
| 12"   | <b>FE27C051202P</b> | pair |
| 18"   | <b>FE27C061802P</b> | pair |

### Adjustable Height Shelf Brackets (in pairs)



| Depth | Part Number         |      |
|-------|---------------------|------|
| 9"    | <b>FE20C040902P</b> | pair |
| 12"   | <b>FE20C051202P</b> | pair |
| 18"   | <b>FE20C061802P</b> | pair |

### Shelf Retaining Lip




| Nominal Length | Painted Steel Part Number | Fits      |      |
|----------------|---------------------------|-----------|------|
| 27"            | <b>FE63M020130N</b>       | 29" Shelf | each |
| 33"            | <b>FE63M020136N</b>       | 35" Shelf | each |
| 42"            | <b>FE63M020142N</b>       | 40" Shelf | each |
| 44"            | <b>FE63M020148N</b>       | 46" Shelf | each |
| 60"            | <b>FE63M020160N</b>       | 58" Shelf | each |
| 72"            | <b>FE63M020172N</b>       | 70" Shelf | each |

| Nominal Length | Stainless Steel Part Number | Fits      |      |
|----------------|-----------------------------|-----------|------|
| 27"            | <b>FE63S020130N</b>         | 29" Shelf | each |
| 33"            | <b>FE63S020136N</b>         | 35" Shelf | each |
| 42"            | <b>FE63S020142N</b>         | 40" Shelf | each |
| 44"            | <b>FE63S020148N</b>         | 46" Shelf | each |
| 60"            | <b>FE63S020160N</b>         | 58" Shelf | each |
| 72"            | <b>FE63S020172N</b>         | 70" Shelf | each |

### LED Task Lighting

**LED Task Lighting** consists of a thin-profile, white, LED strip, 9.5mm high by 33mm deep with a touch dimmer switch. It is UL listed and furnished with 120 VAC to 24 VDC power supply with ten foot power cord. The light includes a magnetic mounting strip for mounting to steel shelving and integral tabs for mounting to wood.



| Part Number        | Length                                    | Fits Under        |
|--------------------|-------------------------------------------|-------------------|
| <b>X-030002-00</b> | 500mm – 19 <sup>5</sup> / <sub>8</sub> "  | 29" - 40" Shelves |
| <b>X-030003-00</b> | 1000mm – 39 <sup>3</sup> / <sub>8</sub> " | 46" - 47" Shelves |

**Note:** Part number is for *single lip only*. Order two lips when lips are required on both front and rear of shelf.


# 4 Free-Standing Table Frames

## Specifications:

**Free-Standing Table Frames** consist of a worksurface frame and a pair of rear adjustable legs. They are available in four lengths, and are adjustable in height from 29" - 36" on 1" increments. The frames are fabricated from powder coated steel with 2" diameter, telescoping legs. Furnished with required installation hardware. Meets SEFA 10, Class 6, Category 4.

**Cabinet Stops** are designed to locate mobile base cabinets in line with the front of the workstation. Fabricated from 16 gauge, powder coated, steel, they are furnished with all required installation hardware.

### Round Free-Standing Table Frames



#### 22" Deep Free-Standing Frame

42" L FE50C362242N  
48" L FE50C362248N  
60" L FE50C362260N  
72" L FE50C362272N

#### 28" Deep Free-Standing Frame


42" L FE50C362842N  
48" L FE50C362848N  
60" L FE50C362860N  
72" L FE50C362872N

#### 34" Deep Free-Standing Frame

42" L FE50C363442N  
48" L FE50C363448N  
60" L FE50C363460N  
72" L FE50C363472N

**ROUND**

### Square Free-Standing Table Frames



#### 22" Deep Free-Standing Frame

42" L FE52C362242N  
48" L FE52C362248N  
60" L FE52C362260N  
72" L FE52C362272N

#### 28" Deep Free-Standing Frame


42" L FE52C362842N  
48" L FE52C362848N  
60" L FE52C362860N  
72" L FE52C362872N

#### 34" Deep Free-Standing Frame

42" L FE52C363442N  
48" L FE52C363448N  
60" L FE52C363460N  
72" L FE52C363472N

**SQUARE**

### Cabinet Stops




#### 28" & 34" Deep Frame Cabinet Stops

for Inset or Overlay Cabinets

42" L FE85C020042N  
48" L FE85C020048N  
60" L FE85C020060N  
72" L FE85C020072N

**Cabinet Stops** can be installed in two positions for the appropriate cabinet positioning. (as pictured)



# Free-Standing Table Worktops

## Specifications:

**Free-Standing Table Worktops** are available in Kemresin or Phenolic Resin. All are 1" thick, have a 1" front overhang, a 1" rear overhang, and are flush with the ends of the worksurface frames. All edges are finished.

These Worktops are used for both Round and Square Tube Free-Standing Tables.

**Powder Coated Steel Back Splash** is 5" high and mounted to the worksurface frame to provide a 4" high splash. Fabricated from 11 gauge powder coated steel, it is furnished complete with installation hardware. Available in black only.

## Free-Standing Table Worktops



### 24" Deep Worktop for use on 22" Free-Standing Table Frames

|       | Kemresin        | Kemresin   | Black Phenolic Resin |
|-------|-----------------|------------|----------------------|
| 42" L | FE96R012442N-__ | BK = Black | FE96T012442N-ABK     |
| 48" L | FE96R012448N-__ | GR = Grey  | FE96T012448N-ABK     |
| 60" L | FE96R012460N-__ | SL = Slate | FE96T012460N-ABK     |
| 72" L | FE96R012472N-__ | PT = Putty | FE96T012472N-ABK     |


### 30" Deep Worktop for use on 28" Free-Standing Table Frames

|       | Kemresin        | Kemresin   | Black Phenolic Resin |
|-------|-----------------|------------|----------------------|
| 42" L | FE96R013042N-__ | BK = Black | FE96T013042N-ABK     |
| 48" L | FE96R013048N-__ | GR = Grey  | FE96T013048N-ABK     |
| 60" L | FE96R013060N-__ | SL = Slate | FE96T013060N-ABK     |
| 72" L | FE96R013072N-__ | PT = Putty | FE96T013072N-ABK     |

### 36" Deep Worktop for use on 34" Free-Standing Table Frames

|       | Kemresin        | Kemresin   | Black Phenolic Resin |
|-------|-----------------|------------|----------------------|
| 42" L | FE96R013642N-__ | BK = Black | FE96T013642N-ABK     |
| 48" L | FE96R013648N-__ | GR = Grey  | FE96T013648N-ABK     |
| 60" L | FE96R013660N-__ | SL = Slate | FE96T013660N-ABK     |
| 72" L | FE96R013672N-__ | PT = Putty | FE96T013672N-ABK     |

## Black Powder Coated Steel Back Splash



|       |              |
|-------|--------------|
| 42" L | FE96M040342N |
| 48" L | FE96M040348N |
| 60" L | FE96M040360N |
| 72" L | FE96M040372N |

# 4 Steel and Wood Cabinets

Storage cabinets can be found in their own catalogs.


## Research Collection Steel Furniture and Signature Series Wood Furniture

### Steel Base Cabinets

- Electrostatically applied, chemical resistant, powder coat finish available in colors to compliment any laboratory environment
- Heavy gauge formed and welded steel cabinet body construction provides strength and durability
- One piece drawer body, radiused all four sides
- Furnished standard with 100 lb. load rated, full-extension, soft-close drawer slides
- One-piece case bottom and horizontal front rail
- Furnished standard with finished back
- Exclusive positive door catch option
- Five-knuckle, semi-concealed, and concealed hinges available


### Wood Base Cabinets

- Multiple coat, chemical resistant finish system formulated to highlight the natural grain patterns
- Dovetailed drawer construction with nine-ply hardwood plywood drawer sides
- 3/4" seven-ply hardwood plywood end panels, bottom and shelf, with solid hardwood facing
- Full extension, ball bearing, 100 lb. load rated, drawer slides, standard on all drawers and pullboards.
- 3/4" seven-ply hardwood plywood finished back dadoed into end panels and bottom
- Institutional type five-knuckle offset door hinges rated for 150 pound loads



# Upper Storage Cabinets

## Upper Cabinet Hanging System



With the use of an upper hanging bracket and a pair of hanging rails, Kewaunee's standard wall cabinets can be mounted on the Enterprise workstation support frame. Two hanging rails are snapped into the support frame slots. The upper hanging bracket bolts to the top of the upper cabinet which then hangs on the upper hanging rail. The lower hanging rail allows the bottom of the cabinet to be securely fastened to the module.


Hanging rails and upper hanging brackets, are finished to match the workstation support frame and are furnished with all required hardware.

Hanging rails are furnished in pairs.

Order the upper hanging bracket in the same length as the upper cabinet. The hanging rails must be the same length as the workstation support frame.


Recommended for use with double-sided support frames, and single-sided support frames installed against a wall.

## Upper Hanging Bracket



| Part Number         | Length |
|---------------------|--------|
| <b>FE68C010430N</b> | 30"    |
| <b>FE68C010436N</b> | 36"    |
| <b>FE68C010442N</b> | 42"    |
| <b>FE68C010448N</b> | 48"    |
| <b>FE68C010460N</b> | 60"    |
| <b>FE68C010472N</b> | 72"    |

## Hanging Rails




| 42" Long Frames       |        | 48" Long Frames     |        |
|-----------------------|--------|---------------------|--------|
| <b>FE69C020142-00</b> | 1 pair | <b>FE69C020148P</b> | 1 pair |
| 60" Long Frames       |        |                     |        |
| <b>FE69C020130L</b>   | 1 pair | Left Mount          |        |
| <b>FE69C020130R</b>   | 1 pair | Right Mount         |        |
| <b>FE69C020160P</b>   | 1 pair | Full Length         |        |
| 72" Long Frames       |        |                     |        |
| <b>FE69C020136L</b>   | 1 pair | Left Mount          |        |
| <b>FE69C020136R</b>   | 1 pair | Right Mount         |        |
| <b>FE69C020172P</b>   | 1 pair | Full Length         |        |

# 4 Ceiling Service Panels


**Features:**

- Designed to integrate into common ceiling grids for delivering plumbing, electrical, and data services to the lab bench
- Fits in standard 2' x 2' and 2' x 4' ceiling grids
- Electrical devices are UL listed.
- Available in all standard powder coat finish colors
- May be customized to meet your project requirements
- Accommodates multiple lab services including:
  - Common gases
  - Air, vacuum, and water
  - High purity and special gases
  - Power and data
- Accepts multiple fitting options including:
  - Standard or fine control valves
  - Pressure gauges and regulators
  - Quick-connect fittings
  - Communication and data outlets
  - Electrical service outlets




# Ceiling Service Panels

## Typical Ceiling Panel Layouts




### Six Station - 2' x 4' panel

Six groups of 4 service cutouts  
Six power receptacle cutouts  
Four data outlet cutouts




### Four Station - 2' x 2' panel

Four groups of 4 service cutouts  
Four power receptacle cutouts  
Four data outlet cutouts




### Four Station - 2' x 4' panel

Four groups of 6 service cutouts  
Four power receptacle cutouts  
Four data outlet cutouts




### Four Station - 2' x 2' panel

Four groups of 3 service cutouts  
Four power receptacle cutouts  
Two data outlet cutouts




### Four Station - 2' x 4' panel

Four groups of 4 service cutouts  
Four power receptacle cutouts  
Four data outlet cutouts




### Two Station - 2' x 2' panel

Two groups of 6 service cutouts  
Four power receptacle cutouts  
Two data outlet cutouts




### Two Station - 2' x 4' panel

Two groups of 6 service cutouts  
Four power receptacle cutouts  
Four data outlet cutouts




### Two Station - 2' x 2' panel

Two groups of 4 service cutouts  
Two power receptacle cutouts  
Two data outlet cutouts



### Two Station - 2' x 4' panel

Two groups of 4 service cutouts  
Two power receptacle cutouts  
Two data outlet cutouts



### One Station - 2' x 2' panel

One group of 6 service cutouts  
One power receptacle cutout  
One data outlet cutout

# 4 Ceiling Service Panel Accessories


## Specifications:

**Ceiling Service Panels** are fabricated from powder coat 16 gauge steel in two sizes to fit standard ceiling grid. They can be punched to accept most electrical, data, or service fittings. (see pages at the end of this catalog for Ceiling Service Panel Technical Data Sheets)

**Fittings and Fixtures must be ordered separately**


**Electrical Fixtures** are furnished complete with device, box and cover plate and are supplied loose, uninstalled.

### 24" x 24" Ceiling Service Panel




| Part Number        | Dimensions (HxDxL) |
|--------------------|--------------------|
| <b>SC10M022424</b> | 1½" x 24" x 24"    |

### 24" x 48" Ceiling Service Panel



| Part Number        | Dimensions (HxDxL) |
|--------------------|--------------------|
| <b>SC10M022448</b> | 1½" x 24" x 48"    |


## Electrical Fixtures



**Twist-Lock  
Single-circuit 20 amp  
0620-1K**

Black


NEMA L5-20 3-wire, specification grade, 120 VAC, 20 amp, single receptacle.



**Data/Comm Outlet  
0723-0B  
0723-0W**

Black  
White

(2) RJ45 Outlets. Fits in 1¼" x 2" cut-out. Cat6




**Dual-circuit  
Twist-Lock 20 amp  
0624-1K**

Black

NEMA L14-20 4-wire, specification grade receptacle. Capable of carrying two 120 VAC, 20 amp circuits.

**Note: Check local electrical codes for compliance.**




**Single-circuit  
Straight Blade, 20 amp  
H581-1K**

Black

NEMA 5-20, 3-wire, hospital grade, polarized duplex receptacle. 120VAC, 20 amp.

**Note: Check local electrical codes for compliance.**

## Plumbing Quick Connect



**Female  
Panel Mount  
¾" IPS Male**


| Part Number      | Key Color |
|------------------|-----------|
| <b>W-9122-K2</b> | Orange    |
| <b>W-9122-K4</b> | Yellow    |
| <b>W-9122-K5</b> | Blue      |
| <b>W-9122-K7</b> | Purple    |
| <b>W-9122-K8</b> | Brown     |



# 2" Configured Technical Data Sheet

## Enterprise Workstation Support Frame **FE3XC8402XXC-...**

Single-sided, 2 inch Post Support Frame without Center Support



### Circuit Type

- ☐ C1 = Single Circuit - Plug exits Left Post
- ☐ C2 = Single Circuit - Plug exits Right Post
- ☐ C3 = Double Circuit - Plug(s) exits Left Post
- ☐ C4 = Double Circuit - Plug(s) exits Right Post

### Data Type

- ☐ D0 = Cover Plates Covering Cutouts
- ☐ D1 = (2) Cat6 Installed Left, Cover Plate Right
- ☐ D2 = (2) Cat6 Installed Right, Cover Plate Left
- ☐ D3 = (2) Cat6 Pre-Wired Left, Cover Plate Right
- ☐ D4 = (2) Cat6 Pre-Wired Right, Cover Plate Left

### Leg Receptacle

- ☐ L0 = Cover Plates Covering Cutouts
- ☐ L1 = Left Post Duplex, Box-Mount, Cover Plate Right
- ☐ L2 = Right Post Duplex, Box-Mount, Cover Plate Left
- ☐ L3 = Both Post Duplex, Box-Mount

### Receptacle Type

- ☐ SK = Spec. Grade - Black
- ☐ SV = Spec. Grade - Ivory
- ☐ SW = Spec. Grade - White
- ☐ SG = Spec. Grade - Grey
- ☐ SR = Spec. Grade - Red
- ☐ HK = Hospital Grade - Black
- ☐ HV = Hospital Grade - Ivory
- ☐ HW = Hospital Grade - White
- ☐ HG = Hospital Grade - Grey
- ☐ HR = Hospital Grade - Red



|                                                                                                              |       |  |          |
|--------------------------------------------------------------------------------------------------------------|-------|--|----------|
| <b>FE3</b> <input type="checkbox"/> <b>C8402</b> <input type="checkbox"/> <input type="checkbox"/> <b>C-</b> | COLOR |  | ITEM NO. |
| PROJECT NAME:                                                                                                |       |  | QUANTITY |


0225

# 2" Configured Technical Data Sheet

5

## Enterprise Workstation Support Frame FE30C8402XXC-...

Single-sided, 2 inch Post Support Frame with Center Divider



### Circuit Type

- ☐ C1 = Single Circuit - Plug exits Left Post
- ☐ C2 = Single Circuit - Plug exits Right Post
- ☐ C3 = Double Circuit - Plug(s) exits Left Post
- ☐ C4 = Double Circuit - Plug(s) exits Right Post

### Data Type

- ☐ D0 = Cover Plates Covering Cutouts
- ☐ D1 = (2) Cat6 Installed Left, Cover Plate Right
- ☐ D2 = (2) Cat6 Installed Right, Cover Plate Left
- ☐ D3 = (2) Cat6 Pre-Wired Left, Cover Plate Right
- ☐ D4 = (2) Cat6 Pre-Wired Right, Cover Plate Left

### Leg Receptacle

- ☐ L0 = Cover Plates Covering Cutouts
- ☐ L1 = Left Post Duplex, Box-Mount, Cover Plate Right
- ☐ L2 = Right Post Duplex, Box-Mount, Cover Plate Left
- ☐ L3 = Both Post Duplex, Box-Mount

### Receptacle Type

- ☐ SK = Spec. Grade - Black
- ☐ SV = Spec. Grade - Ivory
- ☐ SW = Spec. Grade - White
- ☐ SG = Spec. Grade - Grey
- ☐ SR = Spec. Grade - Red
- ☐ HK = Hospital Grade - Black
- ☐ HV = Hospital Grade - Ivory
- ☐ HW = Hospital Grade - White
- ☐ HG = Hospital Grade - Grey
- ☐ HR = Hospital Grade - Red

KEWAUNEE®

|                                                                |       |          |
|----------------------------------------------------------------|-------|----------|
| FE30C8402 <input type="checkbox"/> <input type="checkbox"/> C- | COLOR | ITEM NO. |
| PROJECT NAME:                                                  |       | QUANTITY |

0225

# 5 6" Configured Technical Data Sheet

## Enterprise Workstation Support Frame **FE3XC8406XXC-...** Double-Sided, 6 inch Post Support Frame without Center Support

**Power Cord**

**Services Left Post**  
Pre-piped  
See Catalog for Part Numbers

Indicates Fittings/Fixtures Required  
*Plugged if not selected*

indicate service required

indicate service required

indicate service required

**Electric on Left Post**  
Pre-wired as selected below

120V AC Duplex  
*if selected*

**Left - Data Option - Right**  
Select in Lower Section

**Service Fittings**  
☐ W-0604-00 Ball Valve  
☐ W-0605-00 Needle Valve

**Piping Option (pick one only)**  
☐ Pipe to top of post  
☐ Pipe to top of post with extra:  
☐ 3ft ☐ 4ft ☐ 5ft  
☐ 6ft ☐ 8ft of hose

indicate service required

indicate service required

indicate service required

**Electric on Left Post**  
Pre-wired as selected below

120V AC Duplex  
*if selected*

**Length: (pick one only)**  
☐ 48" ☐ 60" ☐ 72"  
☐ FE30 ☐ FE34

### Circuit Type

- ☐ C5 = Double Circuit - Plug(s) exits Left Post
- ☐ C6 = Double Circuit - Plug(s) exits Right Post
- ☐ C7 = Quad Circuit - Plugs exits Left Post
- ☐ C8 = Quad Circuit - Plugs exits Right Post

### Data Type

- ☐ D0 = Cover Plates Covering Cutouts
- ☐ D5 = (4) Cat6 Installed Left, Cover Plate Right
- ☐ D6 = (4) Cat6 Installed Right, Cover Plate Left
- ☐ D7 = (4) Cat6 Pre-Wired Left, Cover Plate Right
- ☐ D8 = (4) Cat6 Pre-Wired Right, Cover Plate Left

### Leg Receptacle

- ☐ L0 = Cover Plates Covering Cutouts
- ☐ L4 = Left Post Duplex, Flush-Mount, Cover Plate Right
- ☐ L5 = Right Post Duplex, Flush-Mount, Cover Plate Left
- ☐ L6 = Both Post Duplex, Flush-Mount

### Receptacle Type

- ☐ SK = Spec. Grade - Black
- ☐ SV = Spec. Grade - Ivory
- ☐ SW = Spec. Grade - White
- ☐ SG = Spec. Grade - Grey
- ☐ SR = Spec. Grade - Red
- ☐ HK = Hospital Grade - Black
- ☐ HV = Hospital Grade - Ivory
- ☐ HW = Hospital Grade - White
- ☐ HG = Hospital Grade - Grey
- ☐ HR = Hospital Grade - Red



|                                                                                                              |  |       |          |
|--------------------------------------------------------------------------------------------------------------|--|-------|----------|
| <b>FE3</b> <input type="checkbox"/> <b>C8406</b> <input type="checkbox"/> <input type="checkbox"/> <b>C-</b> |  | COLOR | ITEM NO. |
| PROJECT NAME:                                                                                                |  |       | QUANTITY |

0925

# 6" Configured Technical Data Sheet

# 5

## Enterprise Workstation Support Frame FE30C8406XXC-...

Double-Sided, 6 inch Post Support Frame with Center Support

|                                                                                                                                                                                                                                                                                                                                              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |                                                                                                                                                                                                                                                                                                                                              |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>Power Cord</b></p> <p>_____</p>                                                                                                                                                                                                                                                                                                        | <div style="display: flex; justify-content: space-around;"> <div style="border: 1px solid black; padding: 5px; width: 45%;"> <p><b>Left - Data Option</b><br/>Select in Lower Section</p> </div> <div style="border: 1px solid black; padding: 5px; width: 45%;"> <p><b>Data Option - Right</b><br/>Select in Lower Section</p> </div> </div> <div style="display: flex; justify-content: space-between; margin-top: 10px;"> <div style="width: 45%;"> <p><b>Service Fittings</b></p> <p><input type="checkbox"/> W-0604 Ball Valve</p> <p><input type="checkbox"/> W-0605 Needle Valve</p> </div> <div style="width: 45%;"> <p><b>Service Fittings</b></p> <p>W-0604 Ball Valve <input type="checkbox"/></p> <p>W-0605 Needle Valve <input type="checkbox"/></p> </div> </div> <p><b>Piping Option</b> (pick one only)</p> <p><input type="checkbox"/> Pipe to top of post</p> <p><input type="checkbox"/> Pipe to top of post with extra:</p> <p style="margin-left: 20px;"><input type="checkbox"/> 3ft <input type="checkbox"/> 4ft <input type="checkbox"/> 5ft <input type="checkbox"/> 6ft <input type="checkbox"/> 8ft of hose</p> <p><b>Piping/Hose End Option</b> (pick one only)</p> <p><input type="checkbox"/> 3/8" IPS Threaded Female (Standard)</p> <p><input type="checkbox"/> Keyed Male Quick-Connect W-8111-K_</p> | <p><b>Power Cord</b></p> <p>_____</p>                                                                                                                                                                                                                                                                                                        |
| <p><b>Services Left Post</b><br/>Pre-piped to Top of Post<br/>See Catalog for Part Numbers</p> <p>Indicates Fittings/Fixtures Required<br/><i>Plugged if not selected</i></p> <p>_____</p> <p>_____</p> <p>_____</p> <p><b>Electric on Left Post</b><br/>Pre-wired as selected below</p> <p><b>120V AC Duplex</b><br/><i>If selected</i></p> |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | <p><b>Services Left Post</b><br/>Pre-piped to Top of Post<br/>See Catalog for Part Numbers</p> <p>Indicates Fittings/Fixtures Required<br/><i>Plugged if not selected</i></p> <p>_____</p> <p>_____</p> <p>_____</p> <p><b>Electric on Left Post</b><br/>Pre-wired as selected below</p> <p><b>120V AC Duplex</b><br/><i>if selected</i></p> |
| <p><b>Length:</b> (pick one only)</p> <p><input type="checkbox"/> 60" <input type="checkbox"/> 72"</p>                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |                                                                                                                                                                                                                                                                                                                                              |

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>Circuit Type</b></p> <p><input type="checkbox"/> C5 = Double Circuit - Plug(s) exits Left Post</p> <p><input type="checkbox"/> C6 = Double Circuit - Plug(s) exits Right Post</p> <p><input type="checkbox"/> C7 = Quad Circuit - Plugs exit Left Post</p> <p><input type="checkbox"/> C8 = Quad Circuit - Plugs exit Right Post</p> <p><b>Data Type</b></p> <p><input type="checkbox"/> D0 = Cover Plates Covering Cutouts</p> <p><input type="checkbox"/> D5 = (4) Cat6 Installed Left, Cover Plate Right</p> <p><input type="checkbox"/> D6 = (4) Cat6 Installed Right, Cover Plate Left</p> <p><input type="checkbox"/> D7 = (4) Cat6 Pre-Wired Left, Cover Plate Right</p> <p><input type="checkbox"/> D8 = (4) Cat6 Pre-Wired Right, Cover Plate Left</p> <p><b>Leg Receptacle</b></p> <p><input type="checkbox"/> L0 = Cover Plates Covering Cutouts</p> <p><input type="checkbox"/> L4 = Left Post Duplex, Flush-Mount, Cover Plate Right</p> <p><input type="checkbox"/> L5 = Right Post Duplex, Flush-Mount, Cover Plate Left</p> <p><input type="checkbox"/> L6 = Both Post Duplex, Flush-Mount</p> | <p><b>Receptacle Type</b></p> <p><input type="checkbox"/> SK = Spec. Grade - Black</p> <p><input type="checkbox"/> SV = Spec. Grade - Ivory</p> <p><input type="checkbox"/> SW = Spec. Grade - White</p> <p><input type="checkbox"/> SG = Spec. Grade - Grey</p> <p><input type="checkbox"/> SR = Spec. Grade - Red</p> <p><input type="checkbox"/> HK = Hospital Grade - Black</p> <p><input type="checkbox"/> HV = Hospital Grade - Ivory</p> <p><input type="checkbox"/> HW = Hospital Grade - White</p> <p><input type="checkbox"/> HG = Hospital Grade - Grey</p> <p><input type="checkbox"/> HR = Hospital Grade - Red</p> |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|



|                                                                                     |              |                 |                 |
|-------------------------------------------------------------------------------------|--------------|-----------------|-----------------|
| <p><b>FE30C8406</b> <input type="checkbox"/> <input type="checkbox"/> <b>C-</b></p> | <p>COLOR</p> | <p>ITEM NO.</p> | <p>QUANTITY</p> |
| <p>PROJECT NAME:</p>                                                                |              |                 |                 |


0225

Ceiling Service Panel

24 inch x 24 inch

SC10M022424

Plan View of Ceiling Panel



sketch in cutouts required

Order Devices for Panel Cutouts Separately

| Cutout Options:                                                                     |                                         | Size                                                                 | Quantity                        |
|-------------------------------------------------------------------------------------|-----------------------------------------|----------------------------------------------------------------------|---------------------------------|
|  | <b>Electrical Cutout</b><br>0620 & 0624 | <b>Cutout:</b> 2 1/4" x 3 3/8"<br><b>Faceplate:</b> 2 3/4" x 4 1/2"  | <i>sketch in location above</i> |
|  | <b>Data Cutout</b><br>X-020019          | <b>Cutout:</b> 1 1/4" x 2"<br><b>Escutcheon:</b> 1 1/2" x 2 1/4"     | <i>sketch in location above</i> |
|  | <b>Quick-connect Cutout</b><br>W-9122   | <b>Cutout:</b> 1 1/4" Diameter<br><b>Escutcheon:</b> 1 3/4" Diameter | <i>sketch in location above</i> |
|                                                                                     | <b>Custom Cutout 1</b>                  | <b>Cutout Size:</b><br><i>sketch in size above</i>                   | <i>sketch in location above</i> |
|                                                                                     | <b>Custom Cutout 2</b>                  | <b>Cutout Size:</b><br><i>sketch in size above</i>                   | <i>sketch in location above</i> |
|                                                                                     | <b>Custom Cutout 3</b>                  | <b>Cutout Size:</b><br><i>sketch in size above</i>                   | <i>sketch in location above</i> |

|                    |       |  |          |
|--------------------|-------|--|----------|
| <b>SC10M022424</b> | COLOR |  | ITEM NO. |
| PROJECT NAME:      |       |  | QUANTITY |

# Ceiling Panel Technical Data Sheet


5

## Ceiling Service Panel

24 inch x 48 inch

**SC10M022448**

### Plan View of Ceiling Panel



*sketch in cutouts required*

### Order Devices for Panel Cutouts Separately

| Cutout Options:                                                                     |                                         | Size                                                                 | Quantity                        |
|-------------------------------------------------------------------------------------|-----------------------------------------|----------------------------------------------------------------------|---------------------------------|
|  | <b>Electrical Cutout</b><br>0620 & 0624 | <b>Cutout:</b> 2 1/4" x 3 3/8"<br><b>Faceplate:</b> 2 3/4" x 4 1/2"  | <i>sketch in location above</i> |
|  | <b>Data Cutout</b><br>X-020019          | <b>Cutout:</b> 1 1/4" x 2"<br><b>Escutcheon:</b> 1 1/2" x 2 1/4"     | <i>sketch in location above</i> |
|  | <b>Quick-connect Cutout</b><br>W-9122   | <b>Cutout:</b> 1 1/4" Diameter<br><b>Escutcheon:</b> 1 3/4" Diameter | <i>sketch in location above</i> |
|                                                                                     | <b>Custom Cutout 1</b>                  | <b>Cutout Size:</b><br><i>sketch in size above</i>                   | <i>sketch in location above</i> |
|                                                                                     | <b>Custom Cutout 2</b>                  | <b>Cutout Size:</b><br><i>sketch in size above</i>                   | <i>sketch in location above</i> |
|                                                                                     | <b>Custom Cutout 3</b>                  | <b>Cutout Size:</b><br><i>sketch in size above</i>                   | <i>sketch in location above</i> |

**SC10M022448**

PROJECT NAME:

COLOR


ITEM NO.

QUANTITY

Enterprise Workstation Support Frame

Blank, 2 inch Round or Square Post Support Frame

FE3C8402U



FE3C8402U

PROJECT NAME:

COLOR

ITEM NO

QUANTITY


# 2" Blank Frames Technical Design Info

5

## Enterprise Workstation Support Frame

Blank, 2 inch Round or Square Post Support Frame

**FEXXC8402XXU**



SOME EXAMPLES OF MODIFICATIONS REQUIRING A ZZ: (Consult Special Pricing for other modifications)

- Provide hole patterns in each post for edison iii attachment
- Replace standard levelers with stainless levelers
- Add casters (requires counter-balanced fe80 to maintain ul)
- Reposition the standard horizontal raceway
- Add (1) additional horizontal raceway
- Change the single horiz. Raceway to a double-channel horiz. Raceway
- Frame height mod (add or reduce height below lowest service fixture)

- Extend frame length (longer than 72" w/96" max) - centerpost remains
- Reduce frame length (less than 54" w/24" min) - centerpost removed
- Receive a cutout when no cataloged electric or fixture is ordered
- Add (1) leg electric in addition to those offered standard
- Replace leg receptacle with a data device
- Add a non-cataloged electric or data device - **VOIDS UL**
- Stainless construction) - **VOIDS UL**

TO BE DEFINED ON A "U" TECH SHEET:

- Data, elec., fixture c/o positions
- Data, elec. Cord exit
- Number of circuits and positions

- Modified frame dimensions
- Added horiz. Raceway
- Repositioned horiz. Raceway


FOR ANY SPECIAL DESIGN REQUIREMENT THAT IS NOT DESCRIBED ON THIS PAGE, CONTACT YOUR DEALER SUPPORT MANAGER WHO WILL ARRANGE A DESIGN CONSULTATION DISCUSSION

**KEWAUNEE** 0225

Enterprise Workstation Support Frame

Blank, 2 inch Round or Square Post Support Frame  
with Service Chase below Worksurface

FEXXC8402XXL



FE3C8402LL

PROJECT NAME:

COLOR

ITEM NO

QUANTITY


# 2" Blank Frames Technical Design Info

5

## Enterprise Workstation Support Frame

Blank, 2 inch Round or Square Post Support Frame

**FEXXC8402XXL**



SOME EXAMPLES OF MODIFICATIONS REQUIRING A ZZ: (Consult Special Pricing for other modifications)

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>-Provide hole patterns in each post for Edison III attachment</li> <li>-Replace standard levelers with stainless steel levelers</li> <li>-Add casters (requires counter-balanced fe80 to maintain ul)</li> <li>-Add or remove shelf slots for repositioned horizontal service chase</li> <li>-Reposition the standard horizontal raceway</li> <li>-Add (1) additional horizontal raceway</li> <li>-Change the single horiz. Raceway to a double-channel horiz. Raceway</li> <li>-Frame height mod (add or reduce height below lowest service fixture)</li> </ul> | <ul style="list-style-type: none"> <li>-Extend frame length (longer than 72" w/96" max)</li> <li>-Reduce frame length (less than 54" w/24" min)</li> <li>-Receive a cutout when no cataloged electric or fixture is ordered</li> <li>-Add (1) leg electric in addition to those offered standard</li> <li>-Replace leg receptacle with a data device</li> <li>-Add a non-cataloged electric or data device - <b>VOIDS UL</b></li> <li>-Stainless construction) - <b>VOIDS UL</b></li> </ul> |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

TO BE DEFINED ON A "L" TECH SHEET:

|                                                                                                                                                                  |                                                                                                                                                   |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>-Data, elec., fixture c/o positions</li> <li>-Data, elec. Cord exit</li> <li>-Number of circuits and positions</li> </ul> | <ul style="list-style-type: none"> <li>-Modified frame dimensions</li> <li>-Added horiz. Raceway</li> <li>-Repositioned horiz. Raceway</li> </ul> |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|


**KEWAUNEE** 0225

FOR ANY SPECIAL DESIGN REQUIREMENT THAT IS NOT DESCRIBED ON THIS PAGE, CONTACT YOUR DEALER SUPPORT MANAGER WHO WILL ARRANGE A DESIGN CONSULTATION DISCUSSION

Enterprise Workstation Support Frame

Blank, 2 inch Round or Square Post Support Frame  
with Service Chase below Worksurface

FEXXC8402XXD



FE3

C8402

D-

COLOR

ITEM NO

PROJECT NAME:

QUANTITY


# 2" Blank Frames Technical Design Info

5

## Enterprise Workstation Support Frame

Blank, 2 inch Round or Square Post Support Frame

**FEXXC8402XXD**



SOME EXAMPLES OF MODIFICATIONS REQUIRING A ZZ: (Consult Special Pricing for other modifications)

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>-Provide hole patterns in each post for Edison III attachment</li> <li>-Replace standard levelers with stainless steel levelers</li> <li>-Add casters (requires counter-balanced fe80 to maintain ul)</li> <li>-Add or remove shelf slots for repositioned horizontal service chase</li> <li>-Reposition the standard horizontal raceway</li> <li>-Add (1) additional horizontal raceway</li> <li>-Change the single horiz. Raceway to a double-channel horiz. Raceway</li> <li>-Frame height mod (add or reduce height below lowest service fixture)</li> </ul> | <ul style="list-style-type: none"> <li>-Extend frame length (longer than 72" w/96" max)</li> <li>-Reduce frame length (less than 54" w/24" min)</li> <li>-Receive a cutout when no cataloged electric or fixture is ordered</li> <li>-Add (1) leg electric in addition to those offered standard</li> <li>-Replace leg receptacle with a data device</li> <li>-Add a non-cataloged electric or data device - <b>VOIDS UL</b></li> <li>-Stainless construction) - <b>VOIDS UL</b></li> </ul> |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

TO BE DEFINED ON A "D" TECH SHEET:

|                                                                                                                                                                  |                                                                                                                                                   |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>-Data, elec., fixture c/o positions</li> <li>-Data, elec. Cord exit</li> <li>-Number of circuits and positions</li> </ul> | <ul style="list-style-type: none"> <li>-Modified frame dimensions</li> <li>-Added horiz. Raceway</li> <li>-Repositioned horiz. Raceway</li> </ul> |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|


KEWAUNEE 0225

FOR ANY SPECIAL DESIGN REQUIREMENT THAT IS NOT DESCRIBED ON THIS PAGE, CONTACT YOUR DEALER SUPPORT MANAGER WHO WILL ARRANGE A DESIGN CONSULTATION DISCUSSION

Enterprise Workstation Support Frame

Blank, 6 inch Round or Square Post Support Frame

FEXXC8406XXU



FE3

C8406

U

COLOR

ITEM NO

PROJECT NAME:

QUANTITY


# 6" Blank Frames Technical Design Info

5

## Enterprise Workstation Support Frame

Blank, 6 inch Round or Square Post Support Frame

**FEXXC8406XXU**



SOME EXAMPLES OF MODIFICATIONS REQUIRING A ZZ: (Consult Special Pricing for other modifications)

- Provide hole patterns in each post for edison iii attachment
- Replace standard levelers with stainless levelers
- Add casters (requires counter-balanced fe80 to maintain ul)
- Reposition the standard horizontal raceway
- Add (1) additional horizontal raceway
- Change the single horiz. Raceway to a double-channel horiz. Raceway
- Modify frame to receive vacuubrand fixture

- Frame height modification
- Extend frame length (longer than 72" w/96" max) - centerpost remains
- Reduce frame length (less than 54" w/24" min)- centerpost removed
- Receive a cutout when no cataloged electric or fixture is ordered
- Add (1) leg electric in addition to those offered standard
- Replace leg receptacle with a data device
- Add a non-cataloged electric or data device - **VOIDS UL**
- Stainless construction) - **VOIDS UL**

TO BE DEFINED ON A "U" TECH SHEET:

- Data, elec., fixture c/o positions
- Data, elec. Cord exit
- Number of circuits and positions

- Modified frame dimensions
- Added horiz. Raceway
- Repositioned horiz. Raceway


FOR ANY SPECIAL DESIGN REQUIREMENT THAT IS NOT DESCRIBED ON THIS PAGE, CONTACT YOUR DEALER SUPPORT MANAGER WHO WILL ARRANGE A DESIGN CONSULTATION DISCUSSION

**KEWAUNEE** 0225

Enterprise Workstation Support Frame

Blank, 6 inch Round or Square Post Support Frame  
with Service Chase below Worksurface

FEXXC8406XXL



FE3

C8406

L

COLOR

ITEM NO

PROJECT NAME:

QUANTITY


# 6" Blank Frames Technical Design Info

5

## Enterprise Workstation Support Frame

Blank, 6 inch Round or Square Post Support Frame

**FEXXC8406XXL**



SOME EXAMPLES OF MODIFICATIONS REQUIRING A ZZ: (Consult Special Pricing for other modifications)

-Provide hole patterns in each post for Edison III attachment  
-Replace standard levelers with stainless steel levelers  
-Add casters (requires counter-balanced fe80 to maintain ul)  
-Add or remove shelf slots for repositioned horizontal service chase  
-Reposition the standard horizontal raceway  
-Add (1) additional horizontal raceway  
-Change the single horiz. Raceway to a double-channel horiz. Raceway  
-Frame height mod (add or reduce height below lowest service fixture)

-Frame height modification  
-Extend frame length (longer than 72" w/96" max) - centerpost remains  
-Reduce frame length (less than 54" w/24" min)- centerpost removed  
-Receive a cutout when no cataloged electric or fixture is ordered  
-Add (1) leg electric in addition to those offered standard  
-Replace leg receptacle with a data device  
-Add a non-cataloged electric or data device - **VOIDS UL**  
-Stainless construction) - **VOIDS UL**

TO BE DEFINED ON A "L" TECH SHEET:

-Data, elec, fixture c/o positions  
-Data, elec. Cord exit  
-Number of circuits and positions

-Modified frame dimensions  
-Added horiz. Raceway  
-Repositioned horiz. Raceway


FOR ANY SPECIAL DESIGN REQUIREMENT THAT IS NOT DESCRIBED ON THIS PAGE, CONTACT YOUR DEALER SUPPORT MANAGER WHO WILL ARRANGE A DESIGN CONSULTATION DISCUSSION

**KEWAUNEE** 0225

Enterprise Workstation Support Frame

Blank, 6 inch Round or Square Post Support Frame  
with Service Chase below Worksurface

FEXXC8406XXD



FE3

C8406

D

COLOR

ITEM NO

PROJECT NAME:

QUANTITY

0225

EW-04/26-58


# 6" Blank Frames Technical Design Info

5

## Enterprise Workstation Support Frame

Blank, 6 inch Round or Square Post Support Frame

**FEXXC8406XXD**



SOME EXAMPLES OF MODIFICATIONS REQUIRING A ZZ: (Consult Special Pricing for other modifications)

- Provide hole patterns in each post for Edison III attachment
- Replace standard levelers with stainless steel levelers
- Add casters (requires counter-balanced fe80 to maintain ul)
- Add or remove shelf slots for repositioned horizontal service chase
- Reposition the standard horizontal raceway
- Add (1) additional horizontal raceway
- Change the single horiz. Raceway to a double-channel horiz. Raceway
- Frame height mod (add or reduce height below lowest service fixture)

- Frame height modification
- Extend frame length (longer than 72" w/96" max) - centerpost remains
- Reduce frame length (less than 54" w/24" min)- centerpost removed
- Receive a cutout when no cataloged electric or fixture is ordered
- Add (1) leg electric in addition to those offered standard
- Replace leg receptacle with a data device
- Add a non-cataloged electric or data device - **VOIDS UL**
- Stainless construction) - **VOIDS UL**

TO BE DEFINED ON A "D" TECH SHEET:

- Data, elec, fixture c/o positions
- Data, elec. Cord exit
- Number of circuits and positions



- Modified frame dimensions
- Added horiz. Raceway
- Repositioned horiz. Raceway

FOR ANY SPECIAL DESIGN REQUIREMENT THAT IS NOT DESCRIBED ON THIS PAGE, CONTACT YOUR DEALER SUPPORT MANAGER WHO WILL ARRANGE A DESIGN CONSULTATION DISCUSSION

**KEWAUNEE** 0225

# 5 KEWAUNEE® Enterprise Dimensional Data


## 2" Post Single-sided Workstation Support Frame & Worksurface Frame Dimensions



**FE30C8402XX End Elevation**

| Part Number: | A   | B   | C   |
|--------------|-----|-----|-----|
| FE80C3622XXN | 18" | 22" | 24" |
| FE80C3625XXN | 21" | 25" | 27" |
| FE80C3628XXN | 24" | 28" | 30" |
| FE80C3631XXN | 27" | 31" | 33" |
| FE80C3634XXN | 30" | 34" | 36" |


## 6" Post Single and Double-sided Workstation Support Frame & Worksurface Frame Dimensions




| Part Number: | A   | B   | C   | D   |
|--------------|-----|-----|-----|-----|
| FE80C3622XXN | 18" | 20" | 21" | 27" |
| FE80C3625XXN | 21" | 23" | 24" | 30" |
| FE80C3628XXN | 24" | 26" | 27" | 33" |
| FE80C3631XXN | 27" | 29" | 30" | 36" |
| FE80C3634XXN | 30" | 32" | 33" | 39" |

# 5 KEWAUNEE® Enterprise Dimensional Data


## Wall Standard Grounding Details





## 2" Post Top Shelf Dim. (FE28CXXXXXR)



## 6" Post Top Shelf Dim. (FE28CXXXXXR)




## Free-Standing Table Dimension




# Enterprise Worktop to Column Info

5

## 2" Rear Post with Worktops




## 6" Rear Post with Worktops




# Alternate Chase Frame Dim Data

## 2" Post Blank Frame with Lower Service Chase




## 2" Post Blank Frame with Upper and Lower Service Chase




# Alternate Chase Frame Dim Data

5

## 6" Post Blank Frame with Lower Service Chase



## 6" Post Blank Frame with Upper and Lower Service Chase



# ENTERPRISE

## Free-Standing Workstations

Kewaunee Scientific Corporation is dedicated to manufacturing high quality products for the laboratory marketplace. We offer **Total Laboratory Solutions, Innovative Designs, and Technological Expertise** with a worldwide distribution network to provide laboratories that are truly World Class.

### LABORATORY PRODUCTS FROM KEWAUNEE®

ALPHA® – Adaptable Casework System  
ALPHA Columns – Expandable Workstations  
ALPHA Overhead Carrier Systems  
BasikBench® – Simple 4-Leg Benches or Tables  
Element Bench – Efficient Benches with Plumbed Posts  
Enterprise® – Free-Standing Workstations  
Evolution® – Column Based Workstations  
Kemresin® – Epoxy Resin Counters & Sinks  
Kewaunee Matrix® – Educational Stations  
Research Collection® – Steel Casework  
Signature® Series – Wood Casework  
Supreme Air® – Fume Hoods

For our complete product offering, specifications, catalogs, brochures, case studies, and dealer list, go to [kewaunee.com](http://kewaunee.com), or contact your local sales representative directly.



#### Kewaunee Scientific Corporation

P.O. Box 1842 • Statesville, NC 28687-1842  
704-873-7202 • [kscmarketing@kewaunee.com](mailto:kscmarketing@kewaunee.com) • [kewaunee.com](http://kewaunee.com)

#### Kewaunee Labway India Pvt. Ltd.

No. CA-9A, 2nd Floor, Jigani Link Road  
Jigani, Bangalore 562 106, India  
Phone: 9180-27826725 • Fax: 9180-27826724  
Web Site: [kewaunee.in](http://kewaunee.in)

#### Kewaunee Labway Asia Pte. Ltd.

194 Pandan Loop, #06-22  
Pantech Business Hub, Singapore 128 383  
Phone: 65-6773-0288 • Fax: 65-6773-2322



KEWAUNEE®

Kewaunee, Alpha, BasikBench, Enterprise, Evolution, Kemresin, Kewaunee Matrix, Research Collection, Signature, and Supreme Air are registered trademarks of Kewaunee Scientific Corporation