

Cupsinks & Drain Accessories

Cupsinks



SCRO490N-BE Black Epoxy
SCRO490N-GE Grey Epoxy
SCRO490N-PE Putty Epoxy
SCRO490N-SE Slate Epoxy
 Complete with removable strainer.
 6"x 3" inside dimension. 1/2" I.P.S.
 male straight thread outlet.
 (Overall height is 7 1/2")



0491-BP Black Poly
0491-GP Grey Poly
 Complete with removable
 strainer. 6"x 3" inside dimen-
 sion. 1/2" I.P.S. male straight
 thread outlet.
 (Overall height is 8 1/2")



0492-BP Black Poly
0492-GP Grey Poly
 Complete with removable
 strainer. 9"x 3" inside dimen-
 sion. 1/2" I.P.S. male straight
 thread outlet.
 (Overall height is 7")



SCR0493N-BE Black Epoxy
SCR0493N-GE Grey Epoxy
SCR0493N-PE Putty Epoxy
SCR0493N-SE Slate Epoxy
 Complete with removable type
 316 stainless steel wire screen.
 13 3/8" x 4 1/8" x 5 1/4" I.D. 1/2" I.P.S. male
 straight thread outlet.
 (Overall height is 8")

Stainless Steel Cupsinks




0950-0F Stainless Steel
 Made of type 304 stainless
 steel and has integral cross
 bars. 5" I.D. 1/2" I.P.S. male
 straight thread outlet.

Drain Accessories



0464-00 Polyethylene
 Trap and Sink Outlet
 Adjustable riser (tailpiece)
 can be cut to the desired
 length. Sink outlet is complete
 with plug strainer, gasket and
 locknut. Has a 1/2" I.P.S. male
 straight thread outlet. Run-
 ning type trap can be used
 for either "P" or "S" type drain
 connections. 1/2" I.P.S. male
 straight thread outlet.



0469-00 Polyethylene
 Trap and Sink Outlet
 Adjustable riser (tailpiece)
 with 4 1/2" adjustment that
 does not require cutting. The
 sink outlet is complete with
 plug, strainer, gasket and
 locknut. Has a 1/2" I.P.S. male
 straight thread outlet. Trap
 is the bottle "P" type, with
 removable bottom for clea-
 nout. 1/2" I.P.S. male straight
 thread outlet.

Sink Outlet




0482-BP Black Poly
0482-GP Grey Poly
 Complete with locknut. 1/2"
 I.P.S. male straight thread out-
 let. For use with epoxy resin
 tub sinks and troughs.

Sink Overflow



0494BP-04 4" High
0494BP-06 6" High
0494BP-08 8" High
0494BP-10 10" High
 Sink Overflows are designed
 to be used with **0482-BP**
 sink outlets. 1/2" diameter.

St. Steel Sink Drain



0995-00 Stainless Steel
 Crumb Cup Strainer
 Complete with gaskets, lock-
 nut and basket strainer. 1/2"
 male straight thread outlet.

SF-07/26-107

Kemresin Tub Sinks & Drain Troughs

Kemresin Colors

Kemresin® Tub Sinks and Drain Troughs are available in three standard colors: Black, Grey, and Slate.

To specify color, add color code to the end of the part number.

Example: SUR1005N-BK

Color Codes

BK = Black

GR = Grey

PT = Putty

SL = Slate

Kemresin Tub Sinks

Catalog Number	Inside Dimensions (at Top of Sink)			Drain Location		Overall Sink Height
	Length	Width	Height	A	B	
SUR1000N-__	14"	10"	6"	4"	5"	6 ⁷ / ₈ "
SUR1001N-__	16"	8"	7"	4"	4"	8"
SUR1003N-__	18"	15"	11 ³ / ₈ "	4"	4"	12 ⁷ / ₈ "
(1) SUR1003D-__	37"	15"	11"	4 ¹ / ₂ "	4 ¹ / ₂ "	12"
SUR1005N-__	25"	15"	10"	4"	4"	11 ¹ / ₄ "
SUR1006N-__	18"	15"	8"	4"	4"	9"
(2) SUR1006D-__	36 ³ / ₄ "	15"	8 ¹ / ₄ "	4 ¹ / ₂ "	4 ¹ / ₂ "	9"
SUR1007N-__	24"	16"	8"	4"	8"	9 ¹ / ₈ "
SUR1008N-__	16"	12"	8"	4"	6"	9 ¹ / ₈ "
SUR1011N-__	18"	6 ¹ / ₂ "	6"	9"	3 ¹ / ₄ "	7"
(3) SUR1016N-__	23 ³ / ₈ "	13 ³ / ₄ "	13 ³ / ₄ "	11 ¹ / ₈ "	9 ¹ / ₂ "	17 ³ / ₄ "
SUR1021N-__	16"	16"	7 ¹ / ₂ "	4"	4"	8 ³ / ₄ "
SUR1022N-__	21"	18"	10"	4 ¹ / ₂ "	9"	11"
SUR1023N-__	21 ¹ / ₂ "	15 ¹ / ₂ "	7"	4"	4"	8 ³ / ₈ "
SUR1024N-__	21 ¹ / ₂ "	15 ¹ / ₂ "	11"	4"	4"	12 ¹ / ₈ "
(4) SUR1024D-__	44"	15 ¹ / ₂ "	11"	4"	4"	12 ¹ / ₂ "
SUR1026N-__	35 ¹ / ₂ "	19 ¹ / ₂ "	10"	17 ³ / ₄ "	9 ³ / ₄ "	11 ¹ / ₈ "
SUR1027N-__	42"	20"	10"	21"	10"	11 ¹ / ₄ "
SUR1030N-__	16"	12"	8"	8"	6"	9"
SUR1031N-__	24"	16"	8"	12"	8"	9"
SUR1032N-__	29 ³ / ₄ "	14 ³ / ₄ "	6"	14 ⁷ / ₈ "	7 ³ / ₈ "	7"
SUR1034N-__	21 ¹ / ₂ "	19 ¹ / ₄ "	8"	10 ³ / ₄ "	9 ⁵ / ₈ "	9 ¹ / ₈ "
SUR1035N-__	45 ³ / ₄ "	17 ³ / ₈ "	12 ¹ / ₂ "	22 ⁷ / ₈ "	8 ³ / ₄ "	13 ³ / ₈ "
SUR1036N-__	48 ⁵ / ₈ "	24 ¹ / ₂ "	5 ³ / ₈ "	3 ³ / ₄ "	12 ¹ / ₄ "	6"
(5) SUR1038N-__	12" diameter		7 ⁵ / ₈ "	center		8 ³ / ₈ "
(6) SUR1039N-__	18" diameter		8 ¹ / ₈ "	center		8 ⁷ / ₈ "
SURSM01N-__	12"	6"	6"	6"	3"	7"
SURSM02N-__	19 ¹ / ₂ "	16 ¹ / ₂ "	8"	9 ³ / ₄ "	8 ¹ / ₄ "	9"
SURSM03N-__	24"	21"	8"	6"	4 ¹ / ₂ "	9 ¹ / ₈ "
SURSM04N-__	25 ¹ / ₂ "	16 ¹ / ₂ "	19"	4 ¹ / ₂ "	4 ¹ / ₂ "	20 ¹ / ₈ "
SURSM05N-__	28"	18"	10"	4"	4"	11 ¹ / ₈ "
SURSM06N-__	28"	18"	10"	24"	4"	11 ¹ / ₈ "
SURSM07N-__	30"	16"	10"	4"	4"	11 ¹ / ₈ "
SURSM08N-__	30"	16"	10"	26"	4"	11 ¹ / ₈ "

Catalog Number	Inside Dimensions (at Top of Sink)			Drain Location		Overall Sink Height
	Length	Width	Height	A	B	
SURSM09N-__	30"	18"	18"	4 ¹ / ₂ "	4 ³ / ₄ "	18 ³ / ₄ "
SURSM12N-__	8"	8"	8"	4"	4"	8 ⁷ / ₈ "
SURSM13N-__	14"	10"	8"	4"	5"	9"
SURSM14N-__	23 ¹ / ₄ "	14 ¹ / ₄ "	12 ³ / ₄ "	11 ¹ / ₈ "	9 ¹ / ₄ "	13 ⁵ / ₈ "
SURSM15N-__	19 ³ / ₄ "	18"	24"	9 ⁷ / ₈ "	9"	25"
SURSM16N-__	40"	18"	18 ¹ / ₄ "	20"	9"	19 ¹ / ₄ "
SURSM17N-__	29 ³ / ₄ "	23 ¹ / ₂ "	6"	14 ⁷ / ₈ "	11 ³ / ₄ "	7"
SURSM18N-__	18"	15"	15 ¹ / ₄ "	4 ³ / ₄ "	4 ³ / ₄ "	16 ¹ / ₂ "
SURSM19N-__	35 ¹ / ₂ "	19 ¹ / ₂ "	15 ³ / ₈ "	3 ¹ / ₂ "	4 ¹ / ₂ "	16 ¹ / ₂ "
SURSM20N-__	21"	18"	15 ¹ / ₄ "	4 ¹ / ₂ "	4 ¹ / ₂ "	16 ¹ / ₂ "
SURSM21N-__	60 ¹ / ₂ "	22"	7 ¹ / ₄ "	6"	6"	8 ³ / ₈ "
SURSM22N-__	24"	16"	15 ³ / ₈ "	12"	8"	16 ¹ / ₈ "
SURSM23N-__	28"	18"	15 ¹ / ₂ "	5 ¹ / ₂ "	4"	16 ¹ / ₂ "
SURSM25N-__	21"	18"	15 ¹ / ₈ "	10 ¹ / ₂ "	9"	16 ¹ / ₈ "
SURSM26N-__	21"	11 ³ / ₈ "	15 ¹ / ₈ "	10 ¹ / ₂ "	5 ¹ / ₈ "	16 ¹ / ₈ "
SURSM29N-__	36"	24"	12 ¹ / ₄ "	18"	12"	13 ¹ / ₄ "
SURSM31N-__	24"	15 ³ / ₄ "	12 ¹ / ₈ "	4 ³ / ₄ "	8"	13"
SURSM33N-__	24"	18"	11"	4"	9"	11 ⁷ / ₈ "
(7) SURSM34N-__	15"	14"	10 ¹ / ₂ "	7 ¹ / ₂ "	7"	11 ¹ / ₄ "
SURSM36N-__	14"	10"	6"	7"	5"	6 ⁷ / ₈ "
(8) SURSM37N-__	24"	20"	10 ⁵ / ₈ "	6 ³ / ₄ "	10"	11 ¹ / ₄ "
SURSM39N-__	35"	6 ⁵ / ₈ "	4 ³ / ₄ "	17 ¹ / ₂ "	3 ⁵ / ₁₆ "	5 ¹ / ₂ "
SURSM41N-__	18"	15"	8 ¹ / ₄ "	9"	7 ¹ / ₂ "	9"
SURSM42N-__	25"	15"	10 ¹ / ₂ "	12 ¹ / ₂ "	7 ¹ / ₂ "	11 ¹ / ₄ "
SURSM43N-__	21"	18"	10 ³ / ₈ "	10 ¹ / ₂ "	9"	11 ¹ / ₈ "
SURSM44N-__	36"	16"	11"	12"	4"	12"
SURSM45N-__	27 ³ / ₄ "	14 ³ / ₄ "	8 ¹ / ₂ "	10 ¹ / ₂ "	7 ¹ / ₂ "	9 ³ / ₈ "
SURSM47N-__	43"	15 ¹ / ₂ "	11 ¹ / ₈ "	21 ¹ / ₂ "	7 ³ / ₄ "	12"
SURSM49N-__	36"	16"	8 ³ / ₈ "	6"	8"	9"
SURSM50N-__	24"	24"	12 ¹ / ₂ "	12"	12"	13 ¹ / ₄ "
SURSM51N-__	30"	15 ⁷ / ₈ "	16 ³ / ₈ "	6"	7 ⁵ / ₁₆ "	17 ¹ / ₈ "
SURSM52N-__	45 ¹ / ₈ "	11 ¹ / ₂ "	8 ¹ / ₈ "	22 ⁹ / ₁₆ "	5 ³ / ₄ "	9"
SURSM53N-__	31 ¹ / ₂ "	16"	16 ¹ / ₂ "	7 ¹ / ₈ "	8"	17 ¹ / ₄ "
SURSM55N-__	16"	16"	12 ³ / ₈ "	4 ¹ / ₄ "	4 ¹ / ₄ "	13"
SURSM56N-__	42 ¹ / ₄ "	18 ¹ / ₈ "	12 ⁵ / ₈ "	21 ¹ / ₈ "	9 ¹ / ₁₆ "	13 ³ / ₈ "

- (1) Double tub sink. Distance between centerlines of drain outlets is 19³/₈".
- (2) Double tub sink. Distance between centerlines of drain outlets is 19".
- (3) Center table end sink has 24¹/₂" long by 3³/₄" high back splash without trough cutout. Table side of sink is not tapered.
- (4) Double tub sink. Distance between centerlines of drain outlets is 22¹/₂".
- (5) Cylindrical sink -12" inside diameter at top of sink with center drain outlet.
- (6) Cylindrical sink -18" inside diameter at top of sink with center drain outlet.
- (7) Sides have 1¹/₂" taper per foot in height.
- (8) Exterior of one 21" long side is not tapered.

General Notes:

- Sinks with center drain outlets, except **SUR1011N**, can be furnished with a drain recess sized for a crumb cup strainer or garbage disposal when specified.
- Any sink may be furnished with trough cutout. Specify location and size. (Trough location must be at least 1¹/₂" from inside corner and 2" from bottom.)
- Most sinks can be joined together to form multiple tub sinks provided adjoining sides of sinks are the same size. Sinks longer than 30", wider than 24" or deeper than 12" cannot be factory joined, but they can be factory prepared for field assembly. Assemblies of three or more tub sinks should be no longer than 60" overall. For longer assemblies, one factory joint can be made with a second joint prepared for field assembly.

Kemresin Drop-in Sinks



Kemresin Drop-in Sinks are formed in metal molds to provide a lipped one-piece tub with coved corners and bottoms pitched to the drain outlet. They are made from a combination of modified epoxy resins, upgraded for maximum properties with carefully selected curing resins and inert fillers to produce a thermosetting material.

After being removed from the molds, they are cured at high temperatures for several hours to complete the chemical reaction, and ensure maximum chemical, impact and thermal shock resistance. Drain outlets are sized for 1½" sink outlets unless otherwise specified.

They are designed to be mounted in

Kemresin tops in a specially rabbeted cutout so that no support from below is required.

Drop-in Kemresin Sinks




Part Number	I.D.		Depth	O.D.		Drain Location	
	A	B	C	D	E	G	H
SDR1000N-__	13 ⁷ / ₈ "	9 ⁷ / ₈ "	6"	16 ¹ / ₈ "	12 ¹ / ₈ "	7"	5"
SDR1001N-__	8"	16"	7"	10 ¹ / ₄ "	18 ¹ / ₄ "	4"	4"
SDR1003N-__	17 ¹⁵ / ₁₆ "	14 ¹⁵ / ₁₆ "	11"	20 ³ / ₁₆ "	17 ³ / ₁₆ "	3 ¹ / ₂ "	3 ¹ / ₂ "
SDR1005N-__	25"	15"	10"	27 ¹ / ₄ "	17 ¹ / ₄ "	4"	4"
SDR1006N-__	17 ⁷ / ₈ "	14 ⁷ / ₈ "	8"	20 ³ / ₈ "	17 ³ / ₈ "	3 ¹ / ₂ "	3 ¹ / ₂ "
SDR1008N-__	16"	12"	8"	18 ¹ / ₄ "	14 ¹ / ₄ "	4"	6"
SDR1011N-__	18"	6 ¹ / ₂ "	6"	20 ¹ / ₄ "	8 ³ / ₄ "	9"	3 ¹ / ₄ "
SDR1021N-__	16"	16"	7 ¹ / ₂ "	18 ¹ / ₄ "	18 ¹ / ₄ "	4"	4"
(1) SDR1038N-__	12"	diameter	7 ⁵ / ₈ "	14 ¹ / ₄ "	diameter	centered	
(2) SDR1039N-__	18"	diameter	8 ¹ / ₈ "	20 ¹ / ₄ "	diameter	centered	
SDR1040N-__	15 ³ / ₈ "	11 ³ / ₈ "	8"	17 ⁶ / ₈ "	13 ⁶ / ₈ "	4"	4"
SDR1041N-__	23 ³ / ₈ "	15 ³ / ₈ "	8.3"	25 ⁶ / ₈ "	17 ⁶ / ₈ "	4"	4"
SDR1042N-__	28"	15"	11.8"	30 ¹ / ₄ "	17 ¹ / ₄ "	5"	5"
SDR1043N-__	30"	16"	17.8"	32 ¹ / ₄ "	18 ¹ / ₄ "	5"	8"
SDR1044N-__	12"	8"	6"	14 ¹ / ₄ "	10 ¹ / ₄ "	6"	4"
SDR1045N-__	21"	17"	10"	23 ¹ / ₄ "	19 ¹ / ₄ "	4 ¹ / ₂ "	8 ¹ / ₂ "
SDR1046N-__	21 ¹ / ₂ "	15 ¹ / ₂ "	11"	23 ³ / ₄ "	17 ³ / ₄ "	4"	4"
SDR1047N-__	24"	18"	11"	26 ¹ / ₄ "	20 ¹ / ₄ "	4"	4"
(3) SDR1048N-__	17 ¹ / ₂ "	12 ³ / ₈ "	6 ¹ / ₂ "	19 ¹ / ₂ "	14 ⁵ / ₈ "	8 ³ / ₄ "	6 ³ / ₁₆ "
(4) SDR1049N-__	18"	15"	11"	20 ¹ / ₄ "	17 ¹ / ₄ "	4 ¹ / ₂ "	5"
SDR1070N-__	24"	16"	15 ¹ / ₂ "	26 ¹ / ₄ "	18 ¹ / ₄ "	12"	8"

ADA Compliant


SDR1000A-__	13 ⁷ / ₈ "	9 ⁷ / ₈ "	5"	16 ¹ / ₈ "	12 ¹ / ₈ "	4"	4"
SDR1005A-__	25"	15"	4.8"	27 ¹ / ₄ "	17 ¹ / ₄ "	4"	4"
SDR1006A-__	17 ⁷ / ₈ "	14 ⁷ / ₈ "	5"	20 ³ / ₈ "	17 ³ / ₈ "	3 ¹ / ₂ "	3 ¹ / ₂ "

- (1) Cylindrical sink, 12" inside diameter
(2) Cylindrical sink, 18" inside diameter

- (3) Triangular sink used in X92 and X94 assemblies
(4) Stepped sink, 4³/₄" deep at front



PLAN VIEW
SDR1048N-



SECTION VIEW
SDR1049N-

Kemresin Colors

Kemresin® Drop-in Sinks are available in four standard colors: Black, Grey, Putty, and Slate.

To specify color, add color code to the end of the part number.

Example: SDR1005N-BK

Color Codes

BK = Black **PT** = Putty
GR = Grey **SL** = Slate