

CSi Canadian Scientific

The Future of Laboratory Design

QuickShip 

PRODUCT CATALOGUE

IN ASSOCIATION WITH

AMS
AIR MASTER SYSTEMS





Table of Contents

Terms and Conditions.....	Page 3
Colours and Pulls.....	Page 5
Cabinets.....	Page 7
Specialty Cabinets.....	Page 11
Fillers.....	Page 12
Tables and Benches.....	Page 14
Hoods.....	Page 17
Accessories.....	Page 18

Terms and Conditions

Part Numbers

All Quick Ship orders must utilize specified QS part numbers available in the Quick Ship catalog. Configurations outside of what is displayed on the Quick Ship program specification sheets are not available under the Quick Ship lead times.

Purchase Orders

For Quick Ship processing, please submit items on separate purchase order and clearly indicate "AMS Quick Ship". Purchase orders must include correct part numbers, pricing, reference a quote number if applicable, contain accurate contact information, and have a suitable address for material delivery. Do not combine Quick Ship and non Quick Ship items on the same purchase order. If an order is received with both standard and Quick Ship line items, the Quick Ship items will be broken out as a separate sales order and scheduled to ship separately. This additional step may delay processing.

Quantities

Large orders may be subject to limitations on the quantities. Any orders above that quantity will be treated as a standard order and standard lead times will apply.

Pricing

For pricing on Quick Ship items please contact info@canadianscientific.ca.

Lead Times

Quick Ship orders will ship within 20 business days as calculated from the time of sales order confirmation. All ship dates are estimated, and 20 day lead times do not include weekends or holidays.

Change Orders

Change orders and order cancellations cannot be accepted once a Quick Ship order is processed by Canadian Scientific.

Freight and Shipping

Orders will ship complete based on longest lead-time item. Quick Ship orders will ship within 20 business days as calculated from the time of sales order confirmation, and freight will be arranged through the dealer.

Returns

No returns, exchanges, or restocking can be accepted for Quick Ship items. All sales are final.

Terms and Conditions

Disclaimer

Canadian Scientific reserves the right to modify this program at any time.

Warranty

Standard Air Master Systems warranties applies to Quick Ship items.

Powder Coat

Canadian Scientific is currently offering AMS Petal White as our standard Quick Ship powder coat. See below for colour. **Additional colour options coming soon.

Phenolic Resin Surfaces

Surfaces will have finished edges on all visible sides including a chamfer along the top on all sides. All phenolic surfaces are 1" thick and black. See below for colour. Phenolic for shelving brackets is 3/4" thick and black.

Pulls

Standard pull option on Quick Ship is a stainless steel wire pull with magnetic catches. Locks are not an option on Quick Ship Cabinets.




Receive immediate confirmations and shipping quotes within 24 hours. Our reputable network of global manufacturing locations and local fulfillment centers ensures rapid shipment of reliable parts and products.


Colours & Pulls

These reproductions of the standard colors may vary from actual finished products. If color match is critical, actual samples on metal can be supplied upon request. Special match colors can be supplied when specified. Please contact CSi for pricing information.

CASEWORK COLOURS (POWDER COAT)



Petal White



Light Neutral




Model Grey



Wedgewood Blue

PHENOLIC SURFACES COLOUR



Deep Black

PULLS



Stainless Steel Wire Pull




CSi



CASEWORK

Cabinets

Standing height base cabinets are 34.75" high and 23.00" deep. Units are available in various widths as shown. Drawer units feature 4"H toe kicks, stainless steel wire pulls and 100 lb. full extension soft close slides as standard. Door units feature a single shelf as standard. *Tall cabinets feature 5 shelves.



SINGLE DOOR & DRAWER



QS-SB3518-110L
 QS-SB3518-110R
 QS-SB3524-110L
 QS-SB3524-110R

Available Colors: ☐

DOUBLE DOOR & DRAWER



QS-SB3530-220
 QS-SB3536-220
 QS-SB3548-220

Available Colors: ☐

4 FULL DRAWERS



QS-SB3518-040
 QS-SB3524-040

Available Colors: ☐

DOUBLE DOOR



QS-SB3530-200
 QS-SB3536-200
 QS-SB3548-200

Available Colors: ☐

TALL CABINET



QS-SF8436-2SH

Available Colors: ☐

*Includes 5 Shelves

SINGLE DOOR SINK CABINET



QS-SB3212-100-SR*
 QS-SB3512-100-SR
 QS-SB3524-100-SR


NOTE: ONLY AVAILABLE IN RH DOOR

Available Colors: ☐

*ADA HEIGHT (32")

Cabinets

Standing height base cabinets are 34.75" high and 23.00" deep. Units are available in various widths as shown. Drawer units feature 4"H toe kicks, stainless steel wire pulls and 100 lb. full extension soft close slides as standard. Door units feature a single shelf as standard. *Tall cabinets feature 5 shelves.



WALL CABINET SINGLE DOOR



QS-SW3018-1SHL
QS-SW3018-1SHR
QS-SW3024-1SHL
QS-SW3024-1SHR

Available Colors: ☐
*Includes 2 Shelves

DOUBLE DOOR SINK CABINET



QS-SB3530-200-S
QS-SB3536-200-S
QS-SB3548-200-S

Available Colors: ☐

WALL CABINET DOUBLE DOOR



QS-SW3030-2SH
QS-SW3036-2SH

Available Colors: ☐
*Includes 2 Shelves

ADA SINK CABINET



QS-SADAKSC36

Available Colors: ☐

MOBILE DOOR/DRAWER CABINET



QS-SMB2724-110L
QS-SMB2724-110R

Available Colors: ☐

MOBILE DRAWER CABINET



QS-SMB2724-030

Available Colors: ☐

Cabinets

Standing height base cabinets are 34.75" high and 23.00" deep. Units are available in various widths as shown. Drawer units feature 4"H toe kicks, stainless steel wire pulls and 100 lb. full extension soft close slides as standard. Door units feature a single shelf as standard. *Tall cabinets feature 5 shelves.

Quick Ship

Height

Doors

Half Drawers

Width

Drawers

Hinge

QS-SB3518-110L

Quick Ship

Overall Depth

Bracket Depth

Titan Shelf Kit

Width

Bracket Style

QS-WSK-1360-12B

BLIND CORNER CABINETS



QS-SB3548-LHDLPN
QS-SB3548-RHDRPN

Available Colors: ☐

MOBILE CABINET PHENOLIC TOP



QS-PHE-MOD-2423-0.75BC

Available Colors: ☒

END PEDESTAL



QS-SEP3202-23*
QS-SEP3502-23

*ADA Height (32")

Available Colors: ☐

WALL STANDARDS



QS-MTS1136

*Part Number contains 1 Upright

*Slots hold 2 brackets

Available Colors: ☐

SHELF BRACKETS



QS-SE-12

*1 Pair of 12"D Brackets

Available Colors: ☐

WALL SHELVES



QS-WSK-1330-12B
QS-WSK-1336-12B
QS-WSK-1348-12B

Available Colors: ☐

Cabinets

Standing height base cabinets are 34.75" high and 23.00" deep. Units are available in various widths as shown. Drawer units feature 4"H toe kicks, stainless steel wire pulls and 100 lb. full extension soft close slides as standard. Door units feature a single shelf as standard. *Tall cabinets feature 5 shelves.

KNEESPACE RAILS



QS-SKSR0430
QS-SKSR0436

Available Colors: ☐

KNEESPACE PANELS



QS-SKSP3530-2P
QS-SKSP3536-2P
QS-SKSP3530
QS-SKSP3536

Available Colors: ☐



Pictured: 146" L x 60" D Island made with 36" Double Door Base Cabinets x6 (QS-SB3536-200), Phenolic Surface Material, Knee Space rails x2 (QS-SKSR0436), Knee Space Panels x2 (QS-SKSP3536-2P) and C-Fillers x2 (QS-SCF35-12).



Pictured: 24" Mobile Droor/Drawer Cabinet (QS-SMB2724-110R), 60" L x 30" D Titan Table Frame (QS-TF-AL6030-N), and Phenolic Surface Materials.




Pictured: 36" Double Door and Drawer Base Cabinets (QS-SB3536-220) x2, Phenolic Surface Material, 36" Double Door Wall Cabinets x3 (QS-SW3036-2SH), Knee Space Panels (QS-SKSP3536-2P), Knee space Rail (QS-SKSR0436) and L-Fillers x2 (QS-SLF-35-6).



Pictured: 36" Double Door Sink Cabinet (QS-SB3536-200-S), 36" Double Door and Drawer Base Cabinet (QS-SB3536-220), 36" Double Door Wall Cabinet (QS-SW3036-2SH), Phenolic surface material, Bowl Sink (QS-CS-D55), Faucet (QS-FXTDMMIXGNVB-BH-8), and L-Fillers x2 (QS-SLF-35-6).

Specialty Cabinets

Corrosives storage cabinets are reinforced steel with a hi-tech performance polymer alloy coating, which provides 100% coverage against corrosion. Flammable storage cabinets are reinforced steel and insulated. All Flammable cabinets have synchronized self closing doors.



SINGLE DOOR FLAMMABLE



QS-SFSB3524-22SCCL

NOTE: ONLY AVAILABLE IN LH DOOR
*SELF CLOSING DOORS

Available Colors: ☐

SINGLE DOOR ACID CABINET



QS-SASB3524-22R-CORROSIVES

NOTE: ONLY AVAILABLE IN RH DOOR

Available Colors: ☐

DOUBLE DOOR FLAMMABLE



QS-SFSB3530-22SCC
QS-SFSB3536-22SCC
QS-SFSB3548-22SCC

*SELF CLOSING DOORS

Available Colors: ☐

DOUBLE DOOR ACID CABINET





QS-SASB3530-22-CORROSIVES
QS-SASB3536-22-CORROSIVES
QS-SASB3548-22-CORROSIVES

Available Colors: ☐



Fillers

Base cabinet filler panels are available at 35" high. Note: AMS base cabinets are 23.00" deep. (Example: QS-SLF35-6).



WALL CABINET CORNER FILLERS



QS-SWCCF0-0-13

Available Colors: ☐

WALL CABINET FRONT C-FILLERS



QS-SWCFF30-1
QS-SWCFF30-2
QS-SWCFF30-3

Available Colors: ☐

C-FILLERS



QS-SCF35-12
QS-SCF35-1
QS-SCF35-2
QS-SCF35-6

Available Colors: ☐

REAR FILLERS



QS-SLF32-6*
QS-SLF35-6
QS-FSAS-SLF35-6.625

**For ACID and FLAM Cabinets Only.*

**ADA HEIGHT (32") OPTION AVAILABLE*

Available Colors: ☐

FRONT FILLERS



QS-SFFL35-7
QS-SFFL35-1

Available Colors: ☐

CORNER FILLERS



QS-SICF35-5-SJ

**Includes slip joints*

Available Colors: ☐

SLIP JOINTS



QS-SSJ-35

Available Colors: ☐




CSi



TABLES+BENCHES

Tables and Benches

Table Frames (TF) are available 30" deep with standard widths of 60" and 72". As standard, this 1000 lb. capacity table has hanging rails installed for future suspended cabinet installation. Table Frames height adjustments at 1" increments range between 30" to 38".



TITAN TABLE



QS-TF-AL6030-N
QS-TF-AL7230-N

*TOP NOT INCLUDED

Available Colors: ☐

TITAN WORKBENCH



QS-TWB-AL6030-N
QS-TWB-AL7230-N

*TOP NOT INCLUDED

Available Colors: ☐

PHENOLIC SURFACE - TABLE



QS-PHE-TAB-6030-1.00BC
QS-PHE-TAB-7230-1.00BC

Available Colors: ☒

RACEWAY



QS-TC-60RW-5G
QS-TC-72RW-6G

Available Colors: ☐

PHENOLIC SHELVING



QS-PHE-SHE-2813-0.75BC
QS-PHE-SHE-3413-0.75BC
QS-PHE-SHE-4613-0.75BC
QS-PHE-SHE-5615-0.75BC
QS-PHE-SHE-6815-0.75BC

Available Colors: ☒

BENCH SHELVING



QS-TSK-1560-12B
QS-TSK-1572-12B




Available Colors: ☐



Need assistance choosing the right product or part? Give us a call at (226) 780-4793 or email info@canadianscientific.ca.

Tables and Benches

Table Frames (TF) are available 30" deep with standard widths of 60" and 72". As standard, this 1000 lb. capacity table has hanging rails installed for future suspended cabinet installation. Table Frames height adjustments at 1" increments range between 30" to 38".



PHENOLIC SURFACE - BENCH



QS-PHE-TW1-6030-1.00BC
QS-PHE-TW1-7230-1.00BC

Available Colors: ●

Note: Surfaces will have finished edges on all visible sides including a chamfer along the top edge. Phenolic surfaces are 1" thick and black. Phenolic for shelving brackets and mobile cabinet tops is 3/4" thick and black.



Pictured: 60" L x 30" D Titan Workbench (QS-TWB-AL6030-N), 60" L x 30" D Phenolic Surface Material (QS-PHE-TW1-6030-1.00BC) Bench Shelving x2 (QS-TSK-1560-12B)



Pictured: 72" L x 30" D Titan Table (QS-TF-AL7230-N), 72" L x 30" D Phenolic Surface Material (QS-PHE-TW1-7230-1.00BC)

CSi



HOODS+ACCESSORIES

Fume Hoods

The Solution Hood, designed to handle a multitude of applications. Its streamlined posts and radius airfoil, along with a 36-inch height viewing area, offer a softer look than conventional fume hoods. The hood passed NIH, EPA and ASHRAE 110 testing by a third party. The hood can be used in both CV and VAV designs. Epoxy work surface with 2 sink cutout locations is included with GFI receptacles and a light switch as standard.

SOLUTION HOODS



QS-SOL-48
QS-SOL-60
QS-SOL-72

Available Colors:




Pictured: 60" L x 30" D Solution Hood
(QS-SOL-60) with Double Door Base
Cabinets x2 (QS-SB3530-200-S), Air Flow
Alarm (QS-TEL-500), L-Fillers x2 (QS-SLF35-6)



Accessories

Loose fittings and fixtures to round out your project. Sinks, faucets, switches and ceiling services panels to complete or accessorize your space on time and budget. *Note: The CS-D55 and CS-A55 Sinks cannot be placed inside a fume hood.



REMOTE FAUCETS



QS-FXTCWGN
 QS-FXTCWGNVB
 QS-FXTCWDS
 *DS = Down Spout
 *QS-FXTCWGNVB Shown.

VALVES



QS-FXTAIR
 QS-FXTGAS
 QS-FXTNIT
 QS-FXTVAC

ADA DECK MOUNT FAUCETS



QS-FXTDMMIXGNVB-BH-8
 *Holes for Deck Mount fixtures will need to be added in the field.

SWITCHES



QS-BLSW

CEILING SERVICE PANELS



QS-CSPA-24X24-4E8P
 *24x24 Tile.
 (4) Electrical and (8) Plumbing cutouts.
 Available Colors: ☐

HOOD ALARM



QS-TEL500

SINKS - 25"X15" BOWL



QS-CS-D55 - 10" DEEP
 QS-CS-A55 - 5" DEEP
 *A = ADA Depth
 *QS-CS-D55 Shown.
 Available Colors: ☒

CUP SINK - 3"X6"




QS-CS36
 Available Colors: ☒



CSi Canadian Scientific

The Future of Laboratory Design

 info@canadianscientific.ca

 +1 (226) 780-4793

 canadianscientific.ca

AMS

