BASIKBENCH

Simplicity at its Finest






KEWAUNEE

BasikBench®

- Offers the ideal starting point for a wide range of laboratory, cleanroom, and quality control activities
- Delivers a cost effective option for creating a functional workspace
- BasikBench is a value-priced 4-leg height-adjustable work bench





Configure your BasikBench®

Choose your bench, and whether you would like it to arrive knockdown or fully assembled. Then choose your accessories.




FEATURES

- Shelving adjustable in 1" increments
- Allows for 36" mobile storage units
- Allows for up to four service fixtures
- Data ports standard on uprights
- SEFA 10, Category 3 – 1,000 pound capacity
- Top horizontal support channel for fixture hose to rest in
- Bolt together construction
- Ships knock down or assembled
- Worktop height adjustable on 1" increments

BasikBench® 36

Simplicity at its finest, BasikBench 4-leg laboratory bench comes in several sizes with choice of glides or stabilizing casters.



Shown above is an assembled table with top.

Table = ***BB01C363060-CA***
Resin Top = ***BB36R013060-BK***


C = Stabilizing Caster



Minimum Height

27" with caster
25" with glide

G = Glide



Maximum Height

40" with caster
38" with glide

Available in our standard metal paint colors.

12 Black	**463** Iron Ore	**242** Flint Grey	**85** Satin Grey
94 True Blue	**20** Sand Tan	**61** Light Neutral	**78** Snow White
103 Smokey Blue	**90** Amazon	**102** Agave	**104** Sage
87 Mountain Blue	**62** Driftwood	**97** Invigorate	**99** Rapture Red

Disclaimer: On-screen and printed colors are not exact. Use physical samples for accurate selection.

Part Number	Depth	Length
BB01C362436-XX	24	36
BB01C362442-XX	24	42
BB01C362448-XX	24	48
BB01C362454-XX	24	54
BB01C362460-XX	24	60
BB01C362466-XX	24	66
BB01C362472-XX	24	72
BB01C363036-XX	30	36
BB01C363042-XX	30	42
BB01C363048-XX	30	48
BB01C363054-XX	30	54
BB01C363060-XX	30	60
BB01C363066-XX	30	66
BB01C363072-XX	30	72
BB01C363636-XX	36	36
BB01C363642-XX	36	42
BB01C363648-XX	36	48
BB01C363654-XX	36	54
BB01C363660-XX	36	60
BB01C363666-XX	36	66
BB01C363672-XX	36	72

BB01C36 ☐☐☐☐ – ☐☐

Depth Length C/G A/K

C= CASTER
G= GLIDE

A= ASSEMBLED
K= KNOCKDOWN


BasikBench® 84



Shown above is an assembled table with top.

*Table = **BB01C843060-GA***
*Resin Top = **BB84R013060-BK***

BasikBench® 84 expands above the work surface with uprights. Allowing additional accessories and shelves to be mounted above the work surface. Fixtures and data can also be added to each upright. BasikBench® 84 uprights receive only data bezels pre-installed. No data devices are included. Data bezel accepts all Keystone and QuickPort jacks for hard-wired field installations.

C = Stabilizing Caster	G = Glide
	
Minimum Height at Worksurface 27" with caster 25" with glide at Top of Uprights 74.725" with caster 72.725" with glide	Maximum Height at Worksurface 40" with caster 38" with glide at Top of Uprights 87.725" with caster 85.725" with glide

Available in our standard metal paint colors.

12 Black	**463** Iron Ore	**242** Flint Grey	**85** Satin Grey
94 True Blue	**20** Sand Tan	**61** Light Neutral	**78** Snow White
103 Smokey Blue	**90** Amazon	**102** Agave	**104** Sage
87 Mountain Blue	**62** Driftwood	**97** Invigorate	**99** Rapture Red

Disclaimer: On-screen and printed colors are not exact. Use physical samples for accurate selection.

Part Number	Depth	Length
BB01C843036-XX	30	36
BB01C843042-XX	30	42
BB01C843048-XX	30	48
BB01C843054-XX	30	54
BB01C843060-XX	30	60
BB01C843066-XX	30	66
BB01C843072-XX	30	72
BB01C843636-XX	36	36
BB01C843642-XX	36	42
BB01C843648-XX	36	48
BB01C843654-XX	36	54
BB01C843660-XX	36	60
BB01C843666-XX	36	66
BB01C843672-XX	36	72

BB01C84 ☐☐ Depth ☐☐ Length – ☐ C/G ☐ A/K

C= CASTER
G= GLIDE


A= ASSEMBLED
K= KNOCKDOWN

Tops Materials

NOTE: SHIPPED SEPARATELY

KEMRESIN® EPOXY RESIN

Kemresin® is perfectly suited for use in the lab setting. With ultra-smooth non-porous surface appearance, composed of 1 inch thick molded epoxy resin, these tops come in four different colors. Colors include black, grey, putty and slate. With resistance to many acids, alkalies and solvents, this work top is also resistant to both heat and physical use. These Kemresin® tops can be used wherever the highest physical and chemical resistance properties are required.



Grey = GR Slate = SL Putty = PT Black = BK

Disclaimer: On-screen and printed colors are not exact. Use physical samples for accurate selection.

Part Number	Depth	Length
BB36R012436-XX	24	36
BB36R012442-XX	24	42
BB36R012448-XX	24	48
BB36R012454-XX	24	54
BB36R012460-XX	24	60
BB36R012466-XX	24	66
BB36R012472-XX	24	72
BB36R013036-XX	30	36
BB36R013042-XX	30	42
BB36R013048-XX	30	48
BB36R013054-XX	30	54
BB36R013060-XX	30	60
BB36R013066-XX	30	66
BB36R013072-XX	30	72
BB36R013636-XX	36	36
BB36R013642-XX	36	42
BB36R013648-XX	36	48
BB36R013654-XX	36	54
BB36R013660-XX	36	60
BB36R013666-XX	36	66
BB36R013672-XX	36	72

Part Number	Depth	Length
BB84R013036-XX	30	36
BB84R013042-XX	30	42
BB84R013048-XX	30	48
BB84R013054-XX	30	54
BB84R013060-XX	30	60
BB84R013066-XX	30	66
BB84R013072-XX	30	72
BB84R013636-XX	36	36
BB84R013642-XX	36	42
BB84R013648-XX	36	48
BB84R013654-XX	36	54
BB84R013660-XX	36	60
BB84R013666-XX	36	66
BB84R013672-XX	36	72

Depth Length Color

BB36R01 ☐☐☐☐ – ☐☐


Depth Length Color

BB84R01 ☐☐☐☐ – ☐☐

Accessories

NOTE: ACCESSORIES ARE ORDERED AND SHIPPED SEPARATELY

UPPER SHELF



Part Number	Depth	Length
BB21C051236	12	36
BB21C051242	12	42
BB21C051248	12	48
BB21C051254	12	54
BB21C051260	12	60
BB21C051266	12	66
BB21C051227	12	72
BB21C051536	15	36
BB21C051542	15	42
BB21C051548	15	48
BB21C051554	15	54
BB21C051560	15	60
BB21C051566	15	66
BB21C051572	15	72

Shelves are adjustable on 1" increments. The shelf includes a 1-inch rear overhang to align with the work surface overhang.

To order Retaining Rods add **-ROD** to the end of the upper or top shelf part numbers. Rods are ⅜" powder coated steel.

Example: BB21C051566-**ROD**


LOWER SHELF

Part Number	Depth	Length
BB24C010936	9	36
BB24C010942	9	42
BB24C010948	9	48
BB24C010954	9	54
BB24C010960	9	60
BB24C010966	9	66
BB24C010972	9	72
BB24C011236	12	36
BB24C011242	12	42
BB24C011248	12	48
BB24C011254	12	54
BB24C011260	12	60
BB24C011266	12	66
BB24C011272	12	72
BB24C011536	15	36
BB24C011542	15	42
BB24C011548	15	48
BB24C011554	15	54
BB24C011560	15	60
BB24C011566	15	66
BB24C011572	15	72



The lower shelf is installed across the bench's lower support.

TOP SHELF




Part Number	Depth	Length
BB28C040936	9	36
BB28C040942	9	42
BB28C040948	9	48
BB28C040954	9	54
BB28C040960	9	60
BB28C040966	9	66
BB28C040972	9	72
BB28C041236	12	36
BB28C041242	12	42
BB28C041248	12	48
BB28C041254	12	54
BB28C041260	12	60
BB28C041266	12	66
BB28C041272	12	72
BB28C051536	15	36
BB28C051542	15	42
BB28C051548	15	48
BB28C051554	15	54
BB28C051560	15	60
BB28C051566	15	66
BB28C051572	15	72

To order the Retaining Rods shown in the image, add **-ROD** to the end of the upper or top shelf part numbers. Rods are ⅜" powder coated steel.

Example: BB28C051566-**ROD**

SHELF LIPS



Part Number	Length	Part Number
Front Lip		**Rear Lip**
BB63C020136	36	BB63C020236
BB63C020142	42	BB63C020242
BB63C020148	48	BB63C020248
BB63C020154	54	BB63C020254
BB63C020160	60	BB63C020260
BB63C020166	66	BB63C020266
BB63C020172	72	BB63C020272

Accessories

NOTE: ACCESSORIES ARE ORDERED AND SHIPPED SEPARATELY


HANGING RAIL FOR SUSPENDED CABINET

Part Number	Length
BB89C020236	36
BB89C020242	42
BB89C020248	48
BB89C020254	54
BB89C020260	60
BB89C020266	66
BB89C020272	72



This hanging rail allows for the install of suspended cabinets.

CABINET STOP




Part Number	Length
BB85C020236	36
BB85C020242	42
BB85C020248	48
BB85C020254	54
BB85C020260	60
BB85C020266	66
BB85C020272	72

Requires cabinet hanging rail. Not compatible with CPU or Keyboard Holders.

POWER BAR (20 AMP)

Available in Single or Dual Circuit

Part Number	Length	Single or Dual Circuit (S/D)
BB51C030236	36	S
BB51C030242	42	S
BB51C030248	48	S
BB51C030254	54	S
BB51C030260	60	S
BB51C030266	66	S
BB51C030272	72	S
BB52C030236	36	D
BB52C030242	42	D
BB52C030248	48	D
BB52C030254	54	D
BB52C030260	60	D
BB52C030266	66	D
BB52C030272	72	D



Power bars available in dual or single circuit and can be located on the uprights. Dual circuit has two cords. The Single circuit has one cord. Both are 20 amp and come with twist lock plug.

Accessories

NOTE: ACCESSORIES ARE ORDERED AND SHIPPED SEPARATELY

POWER STRIP MOUNTING RAIL

Part Number	Length
BB50C030136	36
BB50C030142	42
BB50C030148	48
BB50C030154	54
BB50C030160	60
BB50C030166	66
BB50C030172	72




Power strip mounting rail allows for user supplied and installed power strip.

LED TASK LIGHTING

BBLED020320 & BBLED020340 Are kits that include a power supply and power supply mounting bracket.

Part Number	Length
BBLED020320	20
BBLED020340	40

LED Task Lighting consists of a thin-profile, white, LED strip, 9.5 mm high by 33mm deep with a touch dimmer switch. It is UL listed and furnished with a 120 VAC to 24 VDC power supply with ten foot power cord. The light includes a magnetic mounting strip for mounting to steel shelving.



Lights listed below do not include power supplies or power supply mounting brackets. These can be purchased as replacement parts.

Part Number	Length	Fits under
X-030002-00	500mm - 19 5/8	29" - 40" Shelves
X-030003-00	1000mm - 39 3/8	46" - 47" Shelves

GLIDES AND CASTERS

Glides and casters come in sets of four.

Part Number
BBCAST2
BBGLIDE

G = Glide






C = Stabilizing Caster

BIN RAIL

Part Number	Length
BB06C030136	36
BB06C030142	42
BB06C030148	48
BB06C030154	54
BB06C030160	60
BB06C030166	66
BB06C030172	72



This bin rail allows for easy installation of bins via a slotted rail. Bins can easily be moved and reorganized.

Note: Bin not included.

Accessories

NOTE: ACCESSORIES ARE ORDERED AND SHIPPED SEPARATELY

MONITOR ARM




This monitor arm rotates 360 degrees at two points. With an adjustable head with universal monitor mount, this arm easily clamps onto any BasikBench.

BB04C360324

LEG ANCHOR

Part Number	Set of
BB05C000001	4



The leg anchor bracket captures the glide on each leg to prevent movement in any direction. The bracket has a mounting hole for an anchor bolt (not included).

KEYBOARD HOLDER




BB02C020230 / BB02C020236

An accessory to hold a standard keyboard and mouse. Select 30 or 36 to match bench depth.


CABLE MANAGEMENT

Part Number	Length
BB07C030136	36
BB07C030142	42
BB07C030148	48
BB07C030154	54
BB07C030160	60
BB07C030166	66
BB07C030172	72



Cable management organizes power and data cables for an organized lab.

CPU HOLDER



BB03C020230 / BB03C020236

The CPU holder attaches via a bracket mount and swivels side to side. The heavy duty holder is adjusted 12" to 20" high and 3.5" to 9" depth.


FIXTURES

Choose up to two fixtures per side. Fixtures can be installed on the assembled or knockdown version, as well as ordered separately. The AIR/VAC/NIT service fixtures include PVC hose at desired length. GAS service fixtures include stainless steel hoses at desired length. Tech sheet required for installed fixtures. Fixture and hose come preassembled terminated in a female swivel.

Hose Type	Part Number	Total Length
AIR/VAC/NIT	BBFXPVHS046-IN	46"
AIR/VAC/NIT	BBFXPVHS096-IN	96"
AIR/VAC/NIT	BBFXPVHS168-IN	168"
GAS	BBFXSSHS046-IN	46"
GAS	BBFXSSHS096-IN	96"
GAS	BBFXSSHS168-IN	168"

-IN only if installed

Length of hose available above the post can be found by subtracting 39 for top and 44 for bottom locations



Storage Cabinets

Enough clearance for our 36" standing height, mobile, cabinets and accessory rail for hanging suspended cabinets. These cabinets are available in many different finishes as found in Kewaunee's Catalogs.





BASIKBENCH

VISIT KEWAUNEE.COM TO FIND YOUR LOCAL DEALER



Kewaunee Scientific Corporation
P.O. Box 1842 ▪ Statesville, NC 28687-1842
704-873-7202 ▪ kscmarketing@kewaunee.com
kewaunee.com

Kewaunee Labway India Pvt. Ltd.
No. CA-9A, 2nd Floor, Jigani Link Road
Jigani, Bangalore 562 106, India
Phone: 9180-27826725 ▪ Fax: 9180-27826724
Web Site: www.kewaunee.in

Kewaunee Labway Asia Pte. Ltd.
194 Pandan Loop, #06-22
Pantech Business Hub, Singapore 128 383
Phone: 65-6773-0288 ▪ Fax: 65-6773-2322

Kewaunee, Alpha, BasikBench, Enterprise, Evolution, Kemresin, Kewaunee Matrix, Research Collection, Signature, and Supreme Air are registered trademarks of Kewaunee Scientific Corporation

KEWAUNEE®

© Copyright 2025

Digital Version 11/14/25